

IBI Molecular Biology agarose products are highly refined and highly purified. Each lot is physically tested and verified for performance in its' designated application. IBI applies strict testing parameters to insure lot to lot consistency. IBI agarose products are offered in a wide variety of package sizes to fit the needs of any research lab. Each bottle of IBI agarose comes packaged with a specification sheet, preparation instructions, and MSDS for your convenience. *If you require a quantity or specific size of agarose that is not listed here, just give us a call at 1-800-253-4942 or send us an e-mail at info@ibisci.com.*

ORDERING INFORMATION

CATALOG#	DESCRIPTION
IB70035	AGAROSE, MOLECULAR BIOLOGY GRADE-25GM
IB70040	AGAROSE, MOLECULAR BIOLOGY GRADE-100GM
IB70041	AGAROSE, MOLECULAR BIOLOGY GRADE-250GM
IB70042	AGAROSE, MOLECULAR BIOLOGY GRADE-500GM
IB70045	AGAROSE, MOLECULAR BIOLOGY GRADE-1KG



PREPARATION OF AN AGAROSE GEL AND ELECTROPHORESIS BUFFER - DNA

- Select the percentage gel necessary to effectively resolve your sample, using the following guideline:

GEL CONCENTRATIONS AND RESOLVING RANGES

Concentration of Agarose in Gel (% w/v)	Efficient Range of Separations of Linear DNA (kb)
0.3%	5-60
0.6%	1-20
0.7%	.8-10
0.9%	.5-7
1.2%	.4-6
1.5%	.2-3
2.0%	.1-2

- Weigh an appropriate quantity of agarose and place it into a 250ml flask.
NOTE: 0.3% equals 0.3gm of agarose per 100ml of gel volume.
- Make 500ml of either 1X TAE or 1X TBE electrophoresis buffer.
- Add ethidium bromide to the diluted electrophoresis buffer to a final concentration of 0.5µg/ml.
- Add 6.6ml of the 1X electrophoresis buffer containing ethidium bromide (made in step 4) per millimeter of gel thickness desired, up to a maximum of 50ml to a flask containing the agarose.
- Make note of the total solution volume, so that a degree of evaporation can be determined and corrected for.
- Heat the agarose slurry in a microwave oven for approximately 90 seconds. Swirl the flask to ensure grains sticking to the walls enter into the solution, heat the solution again for 30-60 seconds if necessary. Swirl flask again to ensure all the agarose has dissolved into the solution.
- Add deionized water to replace any volume lost through evaporation during the heating process. You are ready to pour your gel.

PHYSICAL SPECIFICATIONS

CAS#:	9012-36-6
Moisture:	≤7%
Ash:	≤0.4%
EEO:	.05-13
Sulfate:	≤10%
Clarity @ 1.5% (NTU):	≤3
Gel Strength @ 1.0% (g/cm ²):	≥1200
Gel Strength @ 1.5% (g/cm ²):	≥2500
Gelling Temperature @ 1.5%:	36°C ± 1.5°C
Melting Temperature @ 1.5%:	88°C ± 1.5°C

MOLECULAR BIOLOGY SPECIFICATIONS

DNase assay:	None Detected
RNase assay:	None Detected
Protease assay:	None Detected
Gel analysis assay:	Pass
DNA Resolution:	1000bp (Finely Resolved)
DNA Extraction assay:	Pass
Gel Background:	Very Low

STORAGE

- Store at room temperature
- Protect from moisture
- Light Sensitive

FEATURES

- Extraordinary mechanical resistance insures safe and easy handling
- Pore size can be varied by modifying the gel concentration
- Exceptionally low absorption of staining reagents
- Excellent transparency allows for high visibility
- Greater thermal stability due to high hysteresis

APPLICATIONS

- Nucleic Acid, Analytical and Preparative Electrophoresis
- Blotting
- Protein Electrophoresis, such as Radial Immunodiffusion.

BASE PAIR DATA

1X TAE Buffer	Gel Concentration	1X TBE Buffer
1500-50bp	2.0	1000-50bp
2500-100bp	1.75	2000-100bp
5000-200bp	1.5	3000-150bp
10,000-250bp	1.25	4000-200bp
16,000-300bp	1.0	8000-300bp
20,000-500bp	0.75	12,000-500bp

DID YOU KNOW?

If you are not going to use your gel immediately after preparation, remove it from the casting fixture and place it into a plastic bag or container. Submerge the gel in 1X electrophoresis buffer containing 1mM NaN. Store at +4°C

3:1 Super-Sieve Agarose is designed to have a larger gel network than IBI Ultra Sieve and is recommended for use in the separation of DNA fragments smaller than 1500bp. Gels made with IBI 3:1 Super Sieve have a higher gel strength than competitive products. The gel is exceptionally firm, but still flexible when handled, thus minimizing cracking. IBI 3:1 Super Sieve has the same melting and gelling temperature as our regular molecular biology grade agarose, therefore allowing fast and easy preparation of gels. IBI 3:1 Super Sieve agarose is recommended for all analytical applications, especially when DNA is recovered for subsequent use in enzymatic procedures.

PHYSICAL SPECIFICATIONS

CAS#:	9012-36-6
Moisture:	≤7%
Ash:	≤0.35%
EEO:	≤12
Sulfate:	≤11%
Clarity @ 1.5% (NTU):	≤5
Gel Strength @ 4.0% (g/cm ²):	≥2000-4200
Gelling Temperature @ 1.5%:	≤40.5°C
Melting Temperature @ 1.5%:	≥93.0°C

MOLECULAR BIOLOGY SPECIFICATIONS

DNase assay:	None Detected
RNase assay:	None Detected
Protease assay:	None Detected
Gel analysis assay:	Pass
DNA Extraction assay:	Pass

FEATURES

- DNA Resolution: bands appear sharp and finely resolved
- DNase/RNase activity: None detected
- Gel Background: Very low after EtBr staining
- Blotting: Very good transference for DNA fragments 154-2176bp in 4% gels
- DNA Binding: Very low

STORAGE

- Store at room temperature, Protect from moisture, Light Sensitive

ORDERING INFORMATION

CATALOG# DESCRIPTION

IB70052	AGAROSE, 3:1 SUPER SIEVE-25GM
IB70053	AGAROSE, 3:1 SUPER SIEVE-250GM

RANGES OF SEPARATION

- 2.0% gel 500-1500bp
- 4.0% gel 150-600bp

**These ranges are approximate and have been calculated in 1X TAE buffer



Ultra Sieve Agarose is designed for molecular screening that improves resolution of small DNA fragments and PCR products. IBI Ultra Sieve agarose has excellent gel strength with clear, sharp, finely resolved bands. This makes IBI Ultra Sieve agarose the perfect choice for efficient separation of small DNA fragments and PCR products.

PHYSICAL SPECIFICATIONS

CAS#:	9012-36-6
Moisture:	≤7%
Ash:	≤0.35%
EEO:	≤12
Sulfate:	≤11%
Gel Strength @ 1.5% (g/cm ²):	≥600
Gel Strength @ 3.0% (g/cm ²):	≥1500
Gelling Temperature @ 1.5%:	≤35.5°C
Melting Temperature @ 1.5%:	≥80.0°C

MOLECULAR BIOLOGY SPECIFICATIONS

DNase assay:	None Detected
RNase assay:	None Detected
Protease assay:	None Detected
DNA Extraction assay:	Pass

RANGES OF SEPARATION

- 1.8% gel 400-1200bp
- 3.0% gel 150-800bp
- 4.5% gel 15-400bp

**These ranges are approximate and have been calculated in 1X TAE buffer. To achieve the best resolution of Ultra Sieve gels they should be stored at +4-8°C for 30 minutes prior to use

STORAGE

- Store at room temperature, Protect from moisture, Light Sensitive

ORDERING INFORMATION

CATALOG# DESCRIPTION

IB70054	AGAROSE, ULTRA SIEVE-25GM
IB70055	AGAROSE, ULTRA SIEVE-250GM

FEATURES

- High resolution of short PCR products and DNA fragments
- Improved clarity enhancing visibility
- High gel strength allows use in blotting procedures
- Excellent gel strength, minimizing cracking or breaking during handling



Low Melting Point Agarose is commonly used as a gel matrix when fragment recovery is desired after electrophoresis. IBI low melt agarose is an excellent choice when paired with IBI PCR/Gel Extract Purification Products (see page 70). The properties of low melt agarose promote quick and thorough digestion of gel slices during extraction procedures. IBI low melt agarose is derived by organic synthesis, which generates methoxylate groups from the basic agarose structure. This produces low melting and gelling temperatures when compared to standard agarose. The lower gelling temperature ensures the agarose will be in a liquid state at a temperature range where "in-gel" manipulations can be performed without prior extraction of the DNA or RNA from the gel slice.

ORDERING INFORMATION

CATALOG#	DESCRIPTION
IB70051	AGAROSE, LOW MELTING POINT-25GM
IB70050	AGAROSE, LOW MELTING POINT-50GM
IB70056	AGAROSE, LOW MELTING POINT-100GM
IB70057	AGAROSE, LOW MELTING POINT-250GM
IB70058	AGAROSE, LOW MELTING POINT-500GM
IB70059	AGAROSE, LOW MELTING POINT-1KG

STORAGE

- Store at room temperature, Protect from moisture, Light Sensitive



LOW MELT AGAROSE

PHYSICAL SPECIFICATIONS

CAS#:	9012-36-6
Moisture:	≤7%
Ash:	≤0.4%
EEO:	≤.12
Sulfate:	≤.12%
Clarity @ 1.5% (NTU):	≤4
Gel Strength @ 1.5% (g/cm ²):	≥500
Gelling Temperature @ 1.5%:	24°C - 28°C
Melting Temperature @ 1.5%:	65.5°C

MOLECULAR BIOLOGY SPECIFICATIONS

DNase assay:	None Detected
RNase assay:	None Detected
Protease assay:	None Detected
DNA Extraction assay:	Pass

**LMP agarose will separate & finely resolve molecules >1000bp.

FEATURES

- Lower gel strength than standard agarose, gels are still handled very easily
- Higher clarity (gel transparency) compared to standard agarose
- Greater sieving capacity compared to standard agarose

APPLICATIONS

- Preparative DNA and RNA electrophoresis
- Cloning of tissue culture cells and viral plaque assays

Pulse Field Gel Electrophoresis Agarose is a linear polymer with a high molecular weight, creating gel structures unlike those of standard agarose. This characteristic, along with its very low sulfate content produces a very strong intercatenary interaction, resulting in a gel with a very high gel strength and high exclusion limit.

ORDERING INFORMATION

CATALOG#	DESCRIPTION
IB70060	AGAROSE, PFGE-25GM
IB70061	AGAROSE, PFGE-50GM
IB70062	AGAROSE, PFGE-100GM
IB70063	AGAROSE, PFGE-250GM
IB70064	AGAROSE, PFGE-500GM
IB70065	AGAROSE, PFGE-1KG

STORAGE

- Store at room temperature, Protect from moisture, Light Sensitive



P.F.G.E. AGAROSE

PHYSICAL SPECIFICATIONS

CAS#:	9012-36-6
Moisture:	≤7%
Ash:	≤0.25%
EEO:	≤.12
Sulfate:	≤.12%
Clarity @ 1.5% (NTU):	≤4
Gel Strength @ 1.0% (g/cm ²):	≥1800
Gel Strength @ 1.5% (g/cm ²):	≥3200
Gelling Temperature @ 1.5%:	36°C ± 1.5°C
Melting Temperature @ 1.5%:	88°C ± 1.5°C

MOLECULAR BIOLOGY SPECIFICATIONS

DNase assay:	None Detected
RNase assay:	None Detected
Protease assay:	None Detected
DNA Extraction assay:	Pass

**PFGE agarose will separate & finely resolve molecules >1000bp.

FEATURES

- Extremely high gel strength
- Greater thermal stability when compared to standard agarose
- Exceptionally low absorption of staining agents

APPLICATIONS

- Pulse field gel electrophoresis
- Conventional electrophoresis and blotting
- Cell and enzyme immobilization

TOLL FREE: 1-800-253-4942

IBI InstaPAGE acrylamide solutions are available in 30% and 40% concentrations in the most popular ratios; 19:1, 29:1, and 37.5:1. These ready to use solutions make it fast and easy to prepare polyacrylamide electrophoresis gels. IBI also offers acrylamide and bisacrylamide in individual solution form, which allows the user to mix bisacrylamide solutions in multiple ratios, giving the user greater flexibility while only having to purchase two components.

PHYSICAL SPECIFICATIONS

CAS#:	79-06-1 / 110-26-9
Conductivity:	Max. 10µmho
pH (Acryl: 10%, 0.1M NaCl):	5.5 - 6.5
Purity (Acrylamide):	Min. 98.5%
Purity (Bisacrylamide):	Min. 99.0%

MOLECULAR BIOLOGY SPECIFICATIONS

DNase assay:	None Detected
RNase assay:	None Detected
Protease assay:	None Detected
Gel analysis assay:	Pass
Polymerization assay:	Pass

STORAGE

- Store at +4°C, Keep tightly sealed.

APPLICATIONS

- Polyacrylamide gels for sequencing and protein electrophoresis.

ORDERING INFORMATION



CATALOG#	DESCRIPTION
IB70000	INSTAPAGE ACRYLAMIDE 30% SOLUTION 19:1-500ML
IB70001	INSTAPAGE ACRYLAMIDE 30% SOLUTION 19:1-1L
IB70002	INSTAPAGE ACRYLAMIDE 30% SOLUTION 29:1-500ML
IB70003	INSTAPAGE ACRYLAMIDE 30% SOLUTION 29:1-1L
IB70004	INSTAPAGE ACRYLAMIDE 30% SOLUTION 37.5:1-500ML
IB70005	INSTAPAGE ACRYLAMIDE 30% SOLUTION 37.5:1-1L



PHYSICAL SPECIFICATIONS

CAS#:	79-06-1 / 110-26-9
Conductivity:	Max. 10µmho
pH (Acryl: 10%, 0.1M NaCl):	5.5 - 6.5
Purity (Acrylamide):	Min. 98.5%
Purity (Bisacrylamide):	Min. 99.0%

MOLECULAR BIOLOGY SPECIFICATIONS

DNase assay:	None Detected
RNase assay:	None Detected
Protease assay:	None Detected
Gel analysis assay:	Pass
Polymerization assay:	Pass

STORAGE

- Store at +4°C, Keep tightly sealed.

APPLICATIONS

- Polyacrylamide gels for sequencing and protein electrophoresis. These solutions contain 38% (w/v) acrylamide and 2% (w/v) bisacrylamide for a monomer to crosslink your given ratio. Monomer content (MC) 0.95 and crosslinker 0.05CC.

Use the following formula to determine volume of InstaPAGE (ml) required for your concentration. . .

$$\frac{(MC) \times \text{desired gel concentration} \times \text{desired gel volume}}{40\%}$$

Snowflake	- Frozen Storage Required
Thermometer	- Refrigerated Storage Required
Diamond	- Hazardous Materials

ORDERING INFORMATION



CATALOG#	DESCRIPTION
IB70014	INSTAPAGE ACRYLAMIDE 40% SOLUTION 19:1-500ML
IB70015	INSTAPAGE ACRYLAMIDE 40% SOLUTION 19:1-1L
IB70006	INSTAPAGE ACRYLAMIDE 40% SOLUTION 29:1-500ML
IB70007	INSTAPAGE ACRYLAMIDE 40% SOLUTION 29:1-1L
IB70008	INSTAPAGE ACRYLAMIDE 40% SOLUTION 37.5:1-500ML
IB70009	INSTAPAGE ACRYLAMIDE 40% SOLUTION 37.5:1-1L



Polyacrylamide gels are commonly used in electrophoresis applications, formed by free radical polymerization of acrylamide and bis-acrylamide (crosslinker). IBI offers highly purified acrylamide products in a variety of configurations and sizes. IBI acrylamide and bis-acrylamide powders provide today's researcher with the flexibility of mixing any concentration or ratio of these components to meet their needs. IBI acrylamide and bis-acrylamide products have a purity level greater than 99.9%.

ORDERING INFORMATION

CATALOG#	DESCRIPTION
IB70020	ACRYLAMIDE, SOLID-100GM
IB70024	ACRYLAMIDE, SOLID-500GM
IB70026	ACRYLAMIDE, SOLID-1.5Kg
IB70028	ACRYLAMIDE, SOLID-3.0Kg



PHYSICAL SPECIFICATIONS

CAS#:	79-06-1
Formula Weight:	71.08
Molecular Formula:	C ₃ H ₅ NO
Purity:	Min. 99.9%
Melting Point:	85°C ± 1°C
pH (10%, 0.1M NaCl):	6.0 ± 0.5
Acrylic Acid:	<0.001%
Conductivity:	5µmho
Iron:	<0.0001%
Lead:	<0.0001%

MOLECULAR BIOLOGY SPECIFICATIONS

DNase assay:	None Detected
RNase assay:	None Detected
Protease assay:	None Detected
Gel analysis assay:	Pass
Polymerization assay:	Pass

STORAGE

- Store at room temperature, Protect from moisture, Keep tightly sealed.

APPLICATIONS

- Polyacrylamide gels for sequencing and protein electrophoresis.

ORDERING INFORMATION

CATALOG#	DESCRIPTION
IB70100	BISACRYLAMIDE, SOLID-25GM
IB70102	BISACRYLAMIDE, SOLID-100GM



PHYSICAL SPECIFICATIONS

CAS#:	110-26-9
Formula Weight:	154.17
Molecular Formula:	C ₇ H ₁₀ O ₂ N ₂
Purity:	Min. 99.9%
A ₂₉₀ (1% water)	Max. 0.2
Conductivity (2 water)	<5µmho
Acrylic Acid:	Max. 0.001%
Identification:	Pass
pH (1% 0.1M NaCl) @ 25°C:	Min. 5.0

MOLECULAR BIOLOGY SPECIFICATIONS

DNase assay:	None Detected
RNase assay:	None Detected
Protease assay:	None Detected

STORAGE

- Store at room temperature, Protect from moisture, Keep tightly sealed.

APPLICATIONS

- Polyacrylamide gels for sequencing and protein electrophoresis.

ORDERING INFORMATION

CATALOG#	DESCRIPTION
IB70022	ACRYLAMIDE:BISACRYLAMIDE - 19:1-40GM
IB70023	ACRYLAMIDE:BISACRYLAMIDE - 19:1-200GM
IB70016	ACRYLAMIDE:BISACRYLAMIDE - 29:1-40GM
IB70017	ACRYLAMIDE:BISACRYLAMIDE - 29:1-200GM
IB70018	ACRYLAMIDE:BISACRYLAMIDE - 37.5:1-40GM
IB70019	ACRYLAMIDE:BISACRYLAMIDE - 37.5:1-200GM

PHYSICAL SPECIFICATIONS

See Acrylamide / Bisacrylamide (Above)

MOLECULAR BIOLOGY SPECIFICATIONS

See Acrylamide / Bisacrylamide (Above)

STORAGE

- Store at +4°C, Protect from moisture, Keep tightly sealed, Light sensitive.

APPLICATIONS

- Polyacrylamide gels for sequencing and protein electrophoresis.

Acryliquid-40 and InstaBIS-2 can be used to make 19:1, 29:1, and 37.5:1 bisacrylamide solutions. These products include MSDS, Specification Sheets, and Instructions for mixing different ratios.

ORDERING INFORMATION



CATALOG#	DESCRIPTION
IB70010	ACRYLIQUID 40% SOLUTION-500ML

PHYSICAL SPECIFICATIONS

CAS#:	79-06-1
Formula Weight:	71.08
Molecular Formula:	C ₃ H ₅ NO
Conductivity:	Max. 5µmho
Melting Point:	85°C ± 1°C
pH (@ 25°C):	5.5 - 6.5
Free Acrylic Acid:	<0.005%

MOLECULAR BIOLOGY SPECIFICATIONS

DNase assay:	None Detected
RNase assay:	None Detected
Protease assay:	None Detected
Gel analysis assay:	Pass
Polymerization assay:	Pass

STORAGE

- Store at +4°C, Keep tightly sealed, Protect from moisture.

APPLICATIONS

- Polyacrylamide gels for sequencing and protein electrophoresis. Use reference table below to determine the correct volume of Acryliquid-40 and InstaBIS-2 to make 100ml of X% gel. InstaBIS-2 is a 2% (w/v) solution of ultra-pure bisacrylamide in deionized water.

(ml) Acryliquid-40	(ml) InstaBIS-2
=2.375 (X% gel) 19:1	=2.500 (X% gel) 19:1
=2.417 (X% gel) 29:1	=1.667 (X% gel) 29:1
=2.435 (X% gel) 37.5:1	=1.299 (X% gel) 37.5:1

ORDERING INFORMATION



CATALOG#	DESCRIPTION
IB70012	INSTABIS-2 (2%)-500ML

PHYSICAL SPECIFICATIONS

CAS#:	110-26-9
Formula Weight:	154.17
Molecular Formula:	C ₇ H ₁₀ O ₂ N ₂
Conductivity:	<10µmho
Purity:	>99.0%
A ₂₂₀ (2% water):	<0.4
Acrylic Acid:	<0.001%

MOLECULAR BIOLOGY SPECIFICATIONS

DNase assay:	None Detected
RNase assay:	None Detected
Protease assay:	None Detected
Gel analysis assay:	Pass
Polymerization assay:	Pass

STORAGE

- Store at +4°C, Keep tightly sealed.

APPLICATIONS

- Polyacrylamide gels for sequencing and protein electrophoresis. Use reference table above to determine the correct volume of Acryliquid-40 and InstaBIS-2 to make 100ml of X% gel. InstaBIS-2 is a 2% (w/v) solution of ultra-pure bisacrylamide in deionized water.

ORDERING INFORMATION



CATALOG#	DESCRIPTION
IB70080	AMMONIUM PERSULFATE-100GM

PHYSICAL SPECIFICATIONS

CAS#:	7727-54-0
Formula Weight:	228.19
Molecular Formula:	(NH ₄) ₂ S ₂ O ₈
Purity:	Min. 98.0%
Insolubles:	Max. 0.005%
Residue on Ignition:	Max. 0.05%
Titrateable Free Acid:	Max. 0.04meq/gm
Chloride and Chlorate:	Max. 0.001%
Heavy Metals (Pb):	Max. 0.005%
Iron:	Max. 0.001%
Manganese:	Max. 0.00005%

MOLECULAR BIOLOGY SPECIFICATIONS

Polymerization assay:	Pass
-----------------------	------

STORAGE

- Store at room temperature, Protect from moisture.

APPLICATIONS

- A potent oxidizer that promotes polymerization of acrylamide gels by scavenging dissolved oxygen in the gel solution, thereby accelerating the acrylamide/bisacrylamide reaction.



AMMONIUM PERSULFATE

ORDERING INFORMATION



CATALOG#	DESCRIPTION
IB70120	TEMED-50ML

PHYSICAL SPECIFICATIONS

CAS#:	110-18-9
Formula Weight:	116.21
Molecular Formula:	C ₄ H ₁₀ N ₂
Purity:	Min. 99.0%

MOLECULAR BIOLOGY SPECIFICATIONS

Polymerization assay:	Pass
Gel analysis:	Pass
Boiling Point:	119 - 121°C
Refractive Index:	1.417 - 1.419

STORAGE

- Store at room temperature, Hygroscopic, Keep tightly sealed.

APPLICATIONS

- Used as a catalyst for polymerization of polyacrylamide gels.



TEMED

❄ Snowflake	- Frozen Storage Required
🌡 Thermometer	- Refrigerated Storage Required
⚠ Diamond	- Hazardous Materials

IBI provides a variety of alcohols in several different sizes to fit the needs of the smaller research lab to the high throughput manufacturing facilities. IBI alcohols demonstrate excellent lot-to-lot consistency and purity; whether you require a high-quality ethanol for precipitating nucleic acids, or you need an excellent cleaning reagent, IBI has the product for your application.



ETHANOL FAMILY

ORDERING INFORMATION

CATALOG#	DESCRIPTION
IB15720	ETHANOL (ANHYDROUS ALCOHOL)-500ML
IB15721	ETHANOL (ANHYDROUS ALCOHOL)-1L
IB15724	ETHANOL (ANHYDROUS ALCOHOL)-4L

PHYSICAL SPECIFICATIONS

CAS#:	64-17-5
Formula Weight:	46.07
Molecular Formula:	C ₂ H ₅ N ₄ OH
Ethanol:	95%
Methanol:	5%
Moisture (KF):	Max. 1.0%
Identification (IR):	Pass

MOLECULAR BIOLOGY SPECIFICATIONS

DNase assay:	None Detected
RNase assay:	None Detected

APPLICATIONS

- Ethanol is widely used for precipitating nucleic acids. The nucleic acid precipitate, which is formed in the presence of moderate concentrations of monovalent cations, is recovered by centrifugation and redissolved in an appropriate buffer at the desired concentration.
- 200 Proof Ethanol is denatured with methyl alcohol.

STORAGE

- Store at room temperature inside a flame proof cabinet.



ISOBUTANOL

ORDERING INFORMATION

CATALOG#	DESCRIPTION
IB05120	ISOBUTANOL-500ML

PHYSICAL SPECIFICATIONS

CAS#:	78-83-1
Formula Weight:	74.12
Molecular Formula:	C ₄ H ₁₀ O
Purity:	>99%
A ₂₆₀ :	0.782 - 0.788g/ml
A ₂₈₀ :	0.02%
Residue on evaporation:	<0.005%
Water:	<0.005%

MOLECULAR BIOLOGY SPECIFICATIONS

DNase assay:	None Detected
RNase assay:	None Detected
Ethidium Bromide Extraction assay:	Pass

STORAGE

- Store at room temperature inside a flame proof cabinet.

APPLICATIONS

- Isobutanol is used to extract ethidium bromide from nucleic acid solutions and to concentrate DNA in aqueous solutions. It is also used to overlay polyacrylamide gels to prevent oxidation.

	Snowflake	- Frozen Storage Required
	Thermometer	- Refrigerated Storage Required
	Diamond	- Hazardous Materials

PHYSICAL SPECIFICATIONS

CAS#:	67-56-1
Formula Weight:	32.0
Molecular Formula:	CH ₃ OH
Purity:	Max. 99.8%
ABS (@ 210nm):	Min. 0.80
ABS (@ 220nm):	Min. 0.40
ABS (@ 230nm):	Min. 0.20
ABS (@ 240nm):	Min. 0.10
ABS (@ 260nm):	Min. 0.04
ABS (@ 280-400nm):	Min. 0.01
Carbonyl Compounds:	Max. 1.0%
Color (APHA):	Min. 10
Liquid Chromatography Suitability:	Pass
Residue after evaporation (%):	Min. 0.001
Titrateable acid (meg/gm):	Min. 0.0003
Titrateable base (meg/gm):	Min. 0.0002
Water:	Min. 0.1%



**METHANOL,
HPLC GRADE**

ORDERING INFORMATION

CATALOG#	DESCRIPTION
IB15750	METHANOL, HPLC GRADE-1L

MOLECULAR BIOLOGY SPECIFICATIONS

DNase assay:	None Detected
RNase assay:	None Detected

STORAGE

- Store at room temperature inside a flame proof cabinet.

APPLICATIONS

- Methanol is a high purity, functionally tested reagent suitable for HPLC work. It has low UV absorption and low particulate count to reduce background.

PHYSICAL SPECIFICATIONS

CAS#:	67-56-1
Formula Weight:	32.04
Molecular Formula:	CH ₃ OH
Purity:	Max. 99.8%
Appearance:	Clear, colorless liquid
Specific Gravity:	0.783 - 0.791
Moisture (KF):	Max. 0.1%

STORAGE

- Store at room temperature inside a flame proof cabinet.

MOLECULAR BIOLOGY SPECIFICATIONS

DNase assay:	None Detected
RNase assay:	None Detected

APPLICATIONS

- Ultra-Pure Grade Methanol is used extensively during western transfer procedures. It can also be used for histology and cytology applications.

ORDERING INFORMATION

CATALOG#	DESCRIPTION
IB15755	METHANOL, ULTRA-PURE GRADE-500ML
IB15756	METHANOL, ULTRA-PURE GRADE-1L
IB15757	METHANOL, ULTRA-PURE GRADE-4L



PHYSICAL SPECIFICATIONS

CAS#:	67-63-0
Formula Weight:	60.09
Molecular Formula:	C ₃ H ₇ OH
Purity:	99%
Density:	0.782 - 0.788g/ml
Moisture (KT):	0.02%
Color (APHA):	≤10

MOLECULAR BIOLOGY SPECIFICATIONS

DNase assay:	None Detected
RNase assay:	None Detected

STORAGE

- Store at room temperature inside a flame proof cabinet.

APPLICATIONS

- Isopropanol is widely used for precipitating nucleic acids. The nucleic acid precipitate, which is formed in the presence of moderate concentrations of monovalent cations, is recovered by centrifugation and redissolved in an appropriate buffer at the desired concentrations.

ORDERING INFORMATION

CATALOG#	DESCRIPTION
IB15730	ISOPROPANOL-500ML
IB15735	ISOPROPANOL-1L





In-Vitro culturing of bacterial and animal cells is a routinely performed task in many of today's Life Science Research Laboratories. To ensure the successful growth of target cells it is essential to eliminate non-target bacterial strains and fungi. Antibiotics function in many ways to select against unwanted organisms while maintaining the health and vitality of the desired cells. IBI's antibiotics offer today's researcher excellent potency and lot-to-lot consistency at a price point that will stretch your research dollars!

ORDERING INFORMATION

CATALOG#	DESCRIPTION
IB02010	CARBENICILLIN, DiSODIUM SALT-1GM
IB02020	CARBENICILLIN, DiSODIUM SALT-5GM

STORAGE

- Store at +4°C, Keep tightly sealed, Warm to room temperature prior to use. Store carbenicillin solution in a light-tight container @ -20°C



PHYSICAL SPECIFICATIONS

CAS#:	4800-94-6
Formula Weight:	422.4
Molecular Formula:	C ₁₇ H ₁₆ N ₂ Na ₂ O ₆
Appearance:	White / Off-White Powder
pH (10mg/ml in H ₂ O):	6.0 - 8.0
Solubility:	Clear (50mg/ml in H ₂ O) Solutions are stable for 24 hours at room temperature
Water Content (K.F.):	≤6%

MOLECULAR BIOLOGY SPECIFICATIONS

Antibiotic inhibition/sensitivity assay: Pass

APPLICATIONS

- Carbenicillin is an antibiotic which inhibits bacterial cell wall synthesis activity against gram-positive and gram-negative bacteria, and for cell culture media applications. Recommended use is 100mg/ml

Note: Stock solution of antibiotics dissolved in H₂O should be sterilized by filtration through a 0.22-micron filter.

ORDERING INFORMATION

CATALOG#	DESCRIPTION
IB02030	GENTAMYCIN SULFATE SOLUTION - 20ML TISSUE CULTURE GRADE



STORAGE

- Store at +4°C, Keep tightly sealed.

PHYSICAL SPECIFICATIONS

CAS#:	69-52-3
Potency:	50mg/ml
pH (@ 25°C):	4.5 ± 1.0
Moisture:	Max. 2.0%

MOLECULAR BIOLOGY SPECIFICATIONS

Antibiotic inhibition/sensitivity assay: Pass

APPLICATIONS

- Gentamycin Sulfate Solution is an inhibitor in the growth of gram-positive and gram-negative microorganisms. IBI Gentamycin Sulfate is a broad spectrum antibiotic that is a more effective *in-vitro* bacterial inhibitor than combined streptomycin-penicillin. Gentamycin is non-toxic to tissue culture monolayers, and does not inhibit virus replication. It inhibits the growth of a wide variety of gram-positive and gram-negative microorganisms, including strains resistant to tetracycline, chloramphenicol, kanamycin, and colistin; particularly the strains of *Pseudomonas*, *Proteus*, *Staphylococcus*, and *Streptococcus*. Gentamycin Sulfate inhibits bacterial protein biosynthesis by binding to the 30S subunit of ribosome.

The antibiotics found on pages 34 - 36 are *NOT* for human or veterinary use.

	Snowflake	- Frozen Storage Required
	Thermometer	- Refrigerated Storage Required
	Diamond	- Hazardous Materials

PHYSICAL SPECIFICATIONS

CAS#: 69-52-3
 Formula Weight: 371.39
 Molecular Formula: $C_{16}H_{18}N_3NaO_4S$
 Potency: 845 - 988mcg/mg
 pH (1% water): 9.0 ± 1.0
 Moisture: Max. 2.0%

MOLECULAR BIOLOGY SPECIFICATIONS

Antibiotic inhibition/sensitivity assay: Pass

APPLICATIONS

- Ampicillin inhibits cell wall biosynthesis (peptidoglycan cross-linking). Ampicillin is used at the following working conditions:

	Stock Solution	Working Concentration	
		Stringent Plasmids	Relaxed Plasmids
Ampicillin	50mg/ml in H ₂ O	20µg/ml	60µg/ml

Note: Stock solution of antibiotics dissolved in H₂O should be sterilized by filtration through a 0.22 micron filter.

ORDERING INFORMATION

CATALOG#	DESCRIPTION
IB02040	AMPICILLIN, SODIUM SALT-25GM

STORAGE

- Store at +4°C, Keep tightly sealed, Warm to room temperature prior to use. Store ampicillin solution in a light-tight container @ -20°C



PHYSICAL SPECIFICATIONS

CAS#: 56-75-7
 Formula Weight: 323.13
 Molecular Formula: $C_{11}H_{12}Cl_2N_2O_3$
 Melting Point: 149 - 153°C
 pH (2.5% water): 6.0 ± 1.5
 Purity: 97.0 - 103.0%

MOLECULAR BIOLOGY SPECIFICATIONS

Antibiotic inhibition/sensitivity assay: Pass
 Identification: Pass
 Crystallinity: Pass
 Chromatographic Purity: Pass

APPLICATIONS

- Chloramphenicol is used in the amplification of plasmids and colony hybridizations

	Stock Solution	Working Concentration	
		Stringent Plasmids	Relaxed Plasmids
Chloramphenicol	34mg/ml in ethanol	25µg/ml	170µg/ml

Note: Antibiotics dissolved in ethanol need not be sterilized.

ORDERING INFORMATION

CATALOG#	DESCRIPTION
IB02080	CHLORAMPHENICOL-25GM

STORAGE

- Store at +4°C, Keep tightly sealed, Warm to room temperature prior to use. Store Chloramphenicol solution in a light-tight container @ -20°C



PHYSICAL SPECIFICATIONS

CAS#: 25389-94-0
 Formula Weight: 582.58
 Molecular Formula: $C_{18}H_{36}N_4O_{11}H_2SO_4$
 Potency: >750meg/mg
 pH (1% water): 6.5 - 8.5
 Moisture: Max. 4.0%

MOLECULAR BIOLOGY SPECIFICATIONS

Antibiotic inhibition/sensitivity assay: Pass
 Crystallinity: Pass
 Chromatographic Purity: Pass

Kanamycin Sulfate is derived via microbial fermentation. No animal products or by-products are used in the manufacturing process or the fermentation media.

STORAGE

- Store at +4°C, Keep tightly sealed, Warm to room temperature prior to use. Store Kanamycin Sulfate solution in a light-tight container @ -20°C

ORDERING INFORMATION

CATALOG#	DESCRIPTION
IB02120	KANAMYCIN SULFATE-25GM

APPLICATIONS

- Kanamycin Sulfate is used to bind to 70S ribosomal subunit, inhibits translocation, elicits miscoding bactericidal

	Stock Solution	Working Concentration	
		Stringent Plasmids	Relaxed Plasmids
Kanamycin Sulfate	10mg/ml in H ₂ O	10µg/ml	50µg/ml

Note: Stock solution of antibiotics dissolved in H₂O should be sterilized by filtration through a 0.22 micron filter.

ORDERING INFORMATION

CATALOG# DESCRIPTION

IB02180 STREPTOMYCIN SULFATE-25GM

STORAGE

- Store at +4°C, Keep tightly sealed, Warm to room temperature prior to use. Store Streptomycin Sulfate solution in a light-tight container @ -20°C



PHYSICAL SPECIFICATIONS

CAS#: 3810-74-0
Formula Weight: 1457.38
Molecular Formula: $(C_{21}H_{39}N_7O_{12})_2 \cdot 3H_2SO_4$
Potency: >650mg/mg
pH (1% water): 4.5 - 7.0
Moisture: Max. 4.0%

MOLECULAR BIOLOGY SPECIFICATIONS

Antibiotic inhibition/sensitivity assay: Pass
Crystallinity: Pass
Chromatographic Purity: Pass

APPLICATIONS

- Streptomycin Sulfate is used to bind to S12 protein and 30S ribosomal subunit, inhibits translocation, elicits miscoding bactericidal

	Stock Solution	Working Concentration	
		Stringent Plasmids	Relaxed Plasmids
Streptomycin Sulfate	10mg/ml in H ₂ O	10µg/ml	50µg/ml

Note: Stock solution of antibiotics dissolved in H₂O should be sterilized by filtration through a 0.22 micron filter.

ORDERING INFORMATION

CATALOG# DESCRIPTION

IB02200 TETRACYCLINE HCl-25GM

STORAGE

- Store at +4°C, Keep tightly sealed, Warm to room temperature prior to use. Store Tetracycline HCl solution in a light-tight container @ -20°C



PHYSICAL SPECIFICATIONS

CAS#: 64-75-5
Formula Weight: 480.90
Molecular Formula: $C_{22}H_{24}N_2O_8HCl$
Potency: Min. 900µg/mg
pH (1% water): 2.3 - 0.5
Moisture: Max. 4.0%

MOLECULAR BIOLOGY SPECIFICATIONS

Antibiotic inhibition/sensitivity assay: Pass

APPLICATIONS

- Tetracycline Hydrochloride inhibits protein synthesis (elongation) by preventing the binding of aminoacyl-tRNA to ribosome. Bacteriostatic.

	Stock Solution	Working Concentration	
		Stringent Plasmids	Relaxed Plasmids
Kanamycin Sulfate	5mg/ml in H ₂ O	10µg/ml	50µg/ml

Note: Antibiotics dissolved in ethanol need not be sterilized. Magnesium ions are antagonists of tetracycline. Use media without magnesium salts (e.g. LB medium) for the selection of bacteria resistant to tetracycline.

The antibiotics found on pages 34 - 36 are *NOT* for human or veterinary use.

❄ Snowflake	- Frozen Storage Required
🌡 Thermometer	- Refrigerated Storage Required
💎 Diamond	- Hazardous Materials

BUFFER COMPONENTS

REAGENTS

IBI's buffer components, detergents, and denaturants have been satisfying the needs of Life Science researchers for many years. Excellent purity, convenient package sizes, and competitive pricing makes IBI the perfect choice for your basic Molecular Biology Reagent needs.

ORDERING INFORMATION

CATALOG#	DESCRIPTION
IB15620	AMMONIUM ACETATE-500GM

APPLICATIONS

- Ammonium Acetate is a salt commonly used in the precipitation of nucleic acids. The presence of 2.5M NH₄OAc in the precipitation reaction allows ionically bound contaminants to dissociate from the nucleic acids, allowing pure samples to be recovered. To make 10M Ammonium Acetate Stock Solution, dissolve 385gm of ammonium acetate into 400ml of ddi water. Adjust total volume to 500ml with ddi water. Sterile by filtration.

PHYSICAL SPECIFICATIONS

CAS#:	631-61-8
Molecular Formula:	C ₂ H ₇ NO ₂
Molecular Weight:	77.08
Purity:	Min. 97.0%
Insoluble:	Max. 0.005%
pH (5% water):	6.7 - 7.3
Heavy Metals:	<0.0005%
Iron:	<0.0005%
A ₂₈₀ :	Max. 0.01%
A ₃₅₄ :	Max. 0.01%

MOLECULAR BIOLOGY SPECIFICATIONS

DNase assay:	None Detected
RNase assay:	None Detected
Protease assay:	None Detected

STORAGE

- Store at +4°C, Protect from moisture.

ORDERING INFORMATION

CATALOG#	DESCRIPTION
IB70096	BORIC ACID-2.5Kg



PHYSICAL SPECIFICATIONS

CAS#:	10043-35-3
Molecular Formula:	H ₃ BO ₃
Molecular Weight:	61.83
Purity:	Min. 99.5%
Calcium:	Max. 0.005%
Chloride:	Max. 0.001%
Phosphate:	Max. 0.001%
Sulfate:	Max. 0.01%
Heavy Metals:	Max. 0.001%

MOLECULAR BIOLOGY SPECIFICATIONS

DNase assay:	None Detected
RNase assay:	None Detected
Protease assay:	None Detected

STORAGE

- Store at room temperature, Protect from moisture.

APPLICATIONS

- Boric Acid is an ingredient in Tris-Borate buffer systems, used in high voltage electrophoresis.

ORDERING INFORMATION

CATALOG#	DESCRIPTION
IB70180	EDTA-100GM
IB70182	EDTA-500GM

Buffer	Concentrated Stock Solution (1L)
Tris-acetate (TAE buffer) 50X:	242gm Tris base 57.1ml glacial acetic acid 100ml 0.5M EDTA (pH 8.0)
Tris-borate (TBE buffer) 5X:	54gm Tris base 27.5gm boric acid 20ml 0.5M EDTA (pH 8.0)
Tris-phosphate (TPE buffer) 10X:	108gm Tris base 15.5ml, 85% phosphoric acid (1.679g/ml) 40ml 0.5M EDTA (pH 8.0)

PHYSICAL SPECIFICATIONS

CAS#:	6381-92-6
Molecular Formula:	C ₁₀ H ₁₄ N ₂ O ₈ Na ₂ 2H ₂ O
Molecular Weight:	61.83
Purity:	Min. 99.5%
Calcium:	Max. 0.005%
Chloride:	Max. 0.001%
Phosphate:	Max. 0.001%
Sulfate:	Max. 0.01%
Heavy Metals:	Max. 0.001%

MOLECULAR BIOLOGY SPECIFICATIONS

DNase assay:	None Detected
RNase assay:	None Detected
Protease assay:	None Detected

STORAGE

- Store at room temperature, Protect from moisture.

APPLICATIONS

- EDTA is commonly used in biological and electrophoresis buffer systems. Located in the table to the left are some of the commonly used electrophoresis buffers and their recipes as a concentrated solutions.

BUFFER COMPONENTS

REAGENTS

PHYSICAL SPECIFICATIONS

CAS#:	79-06-1
Molecular Formula:	$C_{10}H_{14}N_2O_8Na \cdot 2H_2O$
Purity:	Min. 99.0%
Insolubles:	Max. 0.01%
Heavy Metals:	<0.005%
A ₂₈₀ (0.5M water):	<0.10%

MOLECULAR BIOLOGY SPECIFICATIONS

DNase assay:	None Detected
RNase assay:	None Detected
Protease assay:	None Detected

STORAGE

- Store at room temperature, Protect from moisture.

APPLICATIONS

- EDTA Solution is commonly used in biological and electrophoresis buffer systems. Located in the table to the right are some of the commonly used electrophoresis buffers and their recipe as a concentrated solution.

ORDERING INFORMATION

CATALOG#	DESCRIPTION
----------	-------------

IB70184	EDTA SOLUTION-100ML
IB70185	EDTA SOLUTION-4x100ML

Buffer		Concentrated Stock Solution (1L)
Tris-acetate (TAE buffer)	50X:	242gm Tris base 57.1ml glacial acetic acid 100ml 0.5M EDTA (pH 8.0)
Tris-borate (TBE buffer)	5X:	54gm Tris base 27.5gm boric acid 20ml 0.5M EDTA (pH 8.0)
Tris-phosphate (TPE buffer)	10X:	108gm Tris base 15.5ml, 85% phosphoric acid (1.679g/ml) 40ml 0.5M EDTA (pH 8.0)

Sambrooks, J., Fritsch, E.F., Maniatis, T (1989) Molecular Cloning A Laboratory Manual Vol. 3 8.23

PHYSICAL SPECIFICATIONS

CAS#:	56-40-6
Molecular Formula:	$C_2H_5NO_2$
Molecular Weight:	75.07
Purity:	Min. 99.5%
A ₂₈₀ (1M water):	<0.10
Heavy Metals:	<0.002%
Solubility (10% water):	Pass
pH (10% water):	5.5 - 6.5
Loss on Drying (%):	≤0.2
Melting Point:	247 ± 2.5°C

MOLECULAR BIOLOGY SPECIFICATIONS

DNase assay:	None Detected
RNase assay:	None Detected
Protease assay:	None Detected

STORAGE

- Store at room temperature, Keep tightly sealed, Protect from moisture.

APPLICATIONS

- Glycine is a neutral amino acid. It is the primary component in polyacrylamide buffers. During the process of separation the glycine forms an electrical front which pulls the macromolecules along.

ORDERING INFORMATION

CATALOG#	DESCRIPTION
----------	-------------

IB70194	GLYCINE-2.5Kg
---------	---------------



PHYSICAL SPECIFICATIONS

CAS#:	7365-45-9
Molecular Formula:	$C_8H_{17}N_2O_4SNa$
Molecular Weight:	260.30
Appearance:	White Powder
Purity:	Min. 99.0
pKa (1X solution @ 25°C):	7.45 - 7.65
Heavy Metals:	<0.0005%

MOLECULAR BIOLOGY SPECIFICATIONS

DNase assay:	None Detected
RNase assay:	None Detected
Protease assay:	None Detected

STORAGE

- Store at room temperature, Protect from moisture.

ORDERING INFORMATION

CATALOG#	DESCRIPTION
----------	-------------

IB01120	HEPES, SODIUM SALT-100GM
---------	--------------------------

APPLICATIONS

- HEPES, Sodium Salt is a zwitterionic buffer ideally suited for cell culture in a range of 10 to 25mM. HEPES has a buffering range from 6.8 to 8.2. HEPES is also a component of Tris-HEPES buffer for polyacrylamide electrophoresis of proteins.

❄ Snowflake	- Frozen Storage Required
🌡 Thermometer	- Refrigerated Storage Required
💎 Diamond	- Hazardous Materials

BUFFER COMPONENTS

REAGENTS

ORDERING INFORMATION

CATALOG#	DESCRIPTION
IB01130	HEPES, FREE ACID-50GM
IB01131	HEPES, FREE ACID-250GM
IB01132	HEPES, FREE ACID-500GM
IB01133	HEPES, FREE ACID-1KG



PHYSICAL SPECIFICATIONS

CAS#:	7365-45-9
Molecular Formula:	$C_8H_{18}N_2O_4S$
Molecular Weight:	238.31
Appearance:	White Powder
Purity:	Min. 99.0
pH (0.5M solution in H_2O):	5.1 - 6.5
Heavy Metals:	Max. 0.001%

MOLECULAR BIOLOGY SPECIFICATIONS

DNase assay:	None Detected
RNase assay:	None Detected
Protease assay:	None Detected

STORAGE

- Store at room temperature, Protect from moisture.

APPLICATIONS

- HEPES, Free Acid is a zwitterionic buffer ideally suited for cell culture in a range of 10 to 25mM. HEPES has a buffering range from 6.8 to 8.2. HEPES is also a component of Tris-HEPES buffer for polyacrylamide electrophoresis of proteins.

ORDERING INFORMATION

CATALOG#	DESCRIPTION
IB70142	TRIS-500GM
IB70144	TRIS-1KG
IB70145	TRIS-5KG

MOLECULAR BIOLOGY SPECIFICATIONS

DNase assay:	None Detected
RNase assay:	None Detected
Protease assay:	None Detected

STORAGE

- Store at room temperature, Protect from moisture.



PHYSICAL SPECIFICATIONS

CAS#:	77-86-1
Molecular Formula:	$C_4H_{11}NO_3$
Molecular Weight:	121.14
Purity:	Min. 99.8%
Moisture:	Max. 1.0%
Insolubles:	Max. 0.005%
A_{280} (1M water):	Max. 0.05
Arsenic:	<0.0005%
Lead:	<0.0001%
Iron:	<0.0001%
Copper:	<0.0001%
Magnesium:	<0.0005%

APPLICATIONS

- Tris is commonly used in biological and electrophoresis buffer systems.

Buffer	Concentrated Stock Solution (1L)
Tris-acetate (TAE buffer) 50X:	242gm Tris base 57.1ml glacial acetic acid 100ml 0.5M EDTA (pH 8.0)
Tris-borate (TBE buffer) 5X:	54gm Tris base 27.5gm boric acid 20ml 0.5M EDTA (pH 8.0)
Tris-phosphate (TPE buffer) 10X:	0.08gm Tris base 15.5ml, 85% phosphoric acid (1.679g/ml) 40ml 0.5M EDTA (pH 8.0)
Tris-glycine (TG buffer) 5X:	15.1gm Tris base 94gm glycine (pH 8.3) 50ml 10% SDS

Sambrooks, J., Fritsch, E.F., Maniatis, T (1989) Molecular Cloning A Laboratory Manual Vol. 3 8.23

	Snowflake	- Frozen Storage Required
	Thermometer	- Refrigerated Storage Required
	Diamond	- Hazardous Materials

BUFFER COMPONENTS

REAGENTS

PHYSICAL SPECIFICATIONS

CAS#:	1185-53-1
Molecular Formula:	$C_4H_{11}NO_3HCl$
Molecular Weight:	157.64
Purity:	Min. 99.5%
Insolubles:	Max. 0.001%
pH (0.1M water @ 25°C):	4.2 - 4.9
Calcium:	<0.0001%
Lead:	<0.0001%
Zinc:	<0.0001%
Copper:	<0.0001%
Magnesium:	<0.0001%
Manganese:	<0.0001%

MOLECULAR BIOLOGY SPECIFICATIONS

DNase assay:	None Detected
RNase assay:	None Detected
Protease assay:	None Detected

STORAGE

- Store at room temperature, Protect from moisture.

APPLICATIONS

- Tris-HCl is commonly used in biological and electrophoresis buffer systems.

ORDERING INFORMATION

CATALOG# DESCRIPTION

IB70162 TRIS-HCl-500GM

Stock Solutions: TE

7.4	Tris-HCl (pH 7.4) 10mM	EDTA (pH 8.0) 1mM
7.6	Tris-HCl (pH 7.6) 10mM	EDTA (pH 8.0) 1mM
8.0	Tris-HCl (pH 8.0) 10mM	EDTA (pH 8.0) 1mM
STE	Tris-HCl (pH 8.0) 10mM	EDTA (pH 8.0 - NaCl) 1mM

PHYSICAL SPECIFICATIONS

CAS#:	151-21-3
Molecular Formula:	$C_{12}H_{25}NaO_4S$
Molecular Weight:	288.38
Purity:	Min. 99.0%
A ₂₆₀ (3% water):	Max. 0.1
A ₂₈₀ (3% water):	Max. 0.1
Chloride:	Max. 0.1%
Phosphate:	Max. 0.0005%
Heavy Metals (Pb):	Max. 0.0005%
Copper:	Max. 0.0005%
% of Straight Chain C12 SDS:	>98

MOLECULAR BIOLOGY SPECIFICATIONS

DNase assay:	None Detected
RNase assay:	None Detected
Protease assay:	None Detected

STORAGE

- Store at room temperature, Protect from moisture.

APPLICATIONS

- SDS is an ionic detergent used to denature proteins in hybridization, nucleic purification, and electrophoresis buffer systems. SDS is also used to dissociate nucleic acid protein complexes in DNA extraction protocols, to disrupt cell membranes and to prepare prehybridization and/or hybridization solutions.
- 10% SDS Solution:** Dissolve 100gm of SDS in 900ml of ddi water. Heat to 68°C to assist dissolution. Adjust the pH to 7.2 by adding a few drops of HCl. Adjust the total volume to 1L with ddi water. Dispense into aliquots.

ORDERING INFORMATION

CATALOG# DESCRIPTION

IB07060 SODIUM DODECYL SULFATE (SDS)-100GM

IB07062 SODIUM DODECYL SULFATE (SDS)-500GM



PHYSICAL SPECIFICATIONS

CAS#:	151-21-3
Molecular Formula:	$C_{12}H_{25}NaO_4S$
Molecular Weight:	288.38
Conductivity:	22800 - 27600 μ mhos
pH (@ 25°C):	5.0 - 8.0

MOLECULAR BIOLOGY SPECIFICATIONS

DNase assay:	None Detected
RNase assay:	None Detected
Protease assay:	None Detected

STORAGE

- Store at room temperature, Protect from moisture.

E-MAIL: INFO@IBISCI.COM

ORDERING INFORMATION

CATALOG# DESCRIPTION

IB07064 20% SDS SOLUTION-100ML

APPLICATIONS

- SDS Solution is an ionic detergent used to denature proteins in hybridization, nucleic purification, and electrophoresis buffer systems. SDS Solution is also used to dissociate nucleic acid protein complexes in DNA extraction protocols, to disrupt cell membranes and to prepare prehybridization and/or hybridization solutions.

BUFFER COMPONENTS

REAGENTS

ORDERING INFORMATION

CATALOG#	DESCRIPTION
IB72060	UREA-500GM
IB72064	UREA-2.5Kg

STORAGE

- Store at room temperature, Protect from moisture.

APPLICATIONS

- Urea is used as a nucleic acid and protein denaturant in polyacrylamide gel electrophoresis. It has a standard concentration of 6 to 8M.

PHYSICAL SPECIFICATIONS

CAS#:	57-13-6
Molecular Formula:	CH ₄ N ₂ O
Molecular Weight:	60.06
Purity:	Min. 99.5%
A ₂₈₀ (8.0M water):	Max. 0.05
Conductivity (5.0M):	Max. 15µmhos
Chloride:	Max. 0.0005%
Iron:	Max. 0.0001%
Heavy Metals:	Max. 0.0005%
Cyanate:	Max. 0.002%

MOLECULAR BIOLOGY SPECIFICATIONS

DNase assay:	None Detected
RNase assay:	None Detected
Protease assay:	None Detected
Sequencing Gel Analysis:	Pass
Biuret Test:	Pass
Solubility (5M water):	Pass

ORDERING INFORMATION

CATALOG#	DESCRIPTION
IB07100	TRITON X-100-100ML



PHYSICAL SPECIFICATIONS

CAS#:	9002-93-1
Molecular Formula:	C ₁₄ H ₂₂ O ₁₂
Molecular Weight:	647.00 (Avg)
pH (5% water):	7.0 ± 1.0
Color APHA:	Max. 100

MOLECULAR BIOLOGY SPECIFICATIONS

DNase assay:	None Detected
RNase assay:	None Detected
Protease assay:	None Detected

STORAGE

- Store at room temperature, Keep tightly sealed.

APPLICATIONS

- Triton X-100 is a non-ionic detergent used to denature cell membranes without denaturing the protein. Triton X-100 breaks up protein aggregates to promote enzymatic activity.

ORDERING INFORMATION

CATALOG#	DESCRIPTION
IB01140	NONIDET P-40 SUBSTITUTE-100ML

APPLICATIONS

- NP-40 is a non-ionic detergent used to denature cell membranes without denaturing the protein. NP-40 detergent is used to solubilize cerebral GABA receptors and solubilization of proteins and lipids.



PHYSICAL SPECIFICATIONS

CAS#:	151-21-3
Appearance:	Colorless, viscous liquid
Density (@ 20°C):	1.05 - 1.07gm/ml
Refractive Index:	1.48 - 1.52

MOLECULAR BIOLOGY SPECIFICATIONS

DNase assay:	None Detected
RNase assay:	None Detected
Protease assay:	None Detected

STORAGE

- Store at room temperature, Protect from moisture. May become turbid after storage, heat to 40°C to clear product.

Snowflake	- Frozen Storage Required
Thermometer	- Refrigerated Storage Required
Diamond	- Hazardous Materials

BUFFER COMPONENTS

REAGENTS

PHYSICAL SPECIFICATIONS

CAS#:	50-01-1
Molecular Formula:	$\text{CH}_5\text{N}_3\text{HCl}$
Molecular Weight:	95.53
Purity:	Max. 0.03
A ₂₆₀ (6M water):	Max. 0.15
A ₂₃₀ (6M water):	4.2 - 4.9
Iron:	Max. 0.0005%
Lead:	Max. 0.0005%

MOLECULAR BIOLOGY SPECIFICATIONS

DNase assay:	None Detected
RNase assay:	None Detected
Protease assay:	None Detected

PHYSICAL SPECIFICATIONS

CAS#:	Blend
Molecular Formula:	$\text{CH}_5\text{N}_3\text{HCl}$
Molecular Weight:	95.53
Solution Components:	Guanidine HCl - 6M EDTA - 25mM
Solution pH:	Adjusted to 7.5

MOLECULAR BIOLOGY SPECIFICATIONS

DNase assay:	None Detected
RNase assay:	None Detected
Protease assay:	None Detected

PHYSICAL SPECIFICATIONS

CAS#:	593-84-0
Molecular Formula:	$\text{CH}_5\text{N}_3\text{CHNS}$
Molecular Weight:	118.16
Purity:	Min. 99.0%
Melting Point:	118 - 121°C
Solubility (70% water):	Clear, haze-free
A ₂₈₀ (70% water):	Max. 0.8
A ₃₀₀ (6M water):	Max. 0.1
A ₄₁₀ (6M water):	Max. 0.1

MOLECULAR BIOLOGY SPECIFICATIONS

DNase assay:	None Detected
RNase assay:	None Detected
Protease assay:	None Detected

STORAGE

- Store at room temperature, Keep tightly sealed, Protect from moisture, Light sensitive.

APPLICATIONS

- Guanidine Thiocyanate is used in the isolation of RNA.

PHYSICAL SPECIFICATIONS

CAS#:	137-16-6
Molecular Formula:	$\text{C}_{18}\text{H}_{35}\text{NNaO}_2$
Molecular Weight:	293.39
Purity:	Min. 95.0%
pH (10% water):	7.5 - 8.5
Sodium Laurate:	Max. 4.0%

MOLECULAR BIOLOGY SPECIFICATIONS

DNase assay:	None Detected
RNase assay:	None Detected
Protease assay:	None Detected

STORAGE

- Store at room temperature, Protect from moisture.

FAX: (563) 690-0490

ORDERING INFORMATION

CATALOG#	DESCRIPTION
IB05080	GUANIDINE HCl-500GM

STORAGE

- Store at room temperature, Protect from moisture.

APPLICATIONS

- Guanidine Hydrochloride is an ionic denaturant that rapidly and effectively denatures most proteins. It is also used in the isolation of RNA.

ORDERING INFORMATION

CATALOG#	DESCRIPTION
IB05085	GUANIDINE HCl SOLUTION-500ML

STORAGE


- Store at room temperature, Protect from moisture.

APPLICATIONS

- Guanidine Hydrochloride Solution is an ionic denaturant that rapidly and effectively denatures most proteins. It is also used in the isolation of RNA.

ORDERING INFORMATION

CATALOG#	DESCRIPTION
IB05100	GUANIDINE THIOCYANATE-500GM



ORDERING INFORMATION

CATALOG#	DESCRIPTION
IB07080	SARKOSYL-100GM

APPLICATIONS

- An ionic detergent used to denature proteins in hybridization, nucleic acid purification, and electrophoresis. It is used in concentrated salt solutions as a detergent because SDS is insoluble in concentrated salt solutions.

BUFFER COMPONENTS

REAGENTS

ORDERING INFORMATION

CATALOG# DESCRIPTION

IB72020	FORMAMIDE, SPECTRAL GRADE-100ML
IB72024	FORMAMIDE, SPECTRAL GRADE-500ML



PHYSICAL SPECIFICATIONS

CAS#:	75-12-7
Molecular Formula:	CH ₃ NO
Molecular Weight:	45.04
Purity:	Min. 99.0%
Conductivity:	Max. 100µmhos
A ₂₈₀ :	Max. 0.08
Iron:	<0.0005%
Copper:	<0.00001%
Lead:	<0.00005%
Zinc:	<0.00005%

MOLECULAR BIOLOGY SPECIFICATIONS

DNase assay:	None Detected
RNase assay:	None Detected
Protease assay:	None Detected

STORAGE

- Store at +4C in the dark, Warm to room temperature prior to opening.
- Formamide is supplied distilled, deionized and packaged under nitrogen. When exposed to oxygen, formamide breaks down into formic acid and ammonia, which can degrade nucleic acid samples. Bottles that have been opened and stored should be deionized prior to use by mixing 1gm of mixed bed resin per 10ml of formamide for 30 minutes.

APPLICATIONS

- Formamide is used as a denaturant in hybridization and electrophoresis applications. Formamide is used to denature DNA and RNA. It separates the two strands of the double helix.

PURITY

- The purity of IBI formamide is determined by way of gas chromatography.

ORDERING INFORMATION

CATALOG# DESCRIPTION

IB72028	FORMAMIDE, BLOT WASH-500ML
---------	----------------------------



PHYSICAL SPECIFICATIONS

CAS#:	75-12-7
Molecular Formula:	CH ₃ NO
Molecular Weight:	45.04
Purity:	Min. 99.0%
pH:	7.1 ± 0.5
Conductivity:	Max. 350µmhos
A ₂₆₀ :	<2.5
A ₂₈₀ :	<0.5
Density:	130g/ml

MOLECULAR BIOLOGY SPECIFICATIONS

DNase assay:	None Detected
RNase assay:	None Detected
Protease assay:	None Detected

STORAGE

- Store at +4C in the dark, Warm to room temperature prior to opening.
- Formamide is supplied distilled, deionized and packaged under nitrogen. When exposed to oxygen, formamide breaks down into formic acid and ammonia, which can degrade nucleic acid samples. Bottles that have been opened and stored should be deionized prior to use by mixing 1gm of mixed bed resin per 10ml of formamide for 30 minutes.

APPLICATIONS

- Blot Washing Formamide is suitable for blot washing procedures to remove background hybridization.

❄ Snowflake	- Frozen Storage Required
🌡 Thermometer	- Refrigerated Storage Required
💎 Diamond	- Hazardous Materials

DYES & STAINS

REAGENTS

IBI's ready-to-use 6X loading dye is ideally suited for agarose electrophoresis of DNA, RNA, or nucleic acids. Its ultra pure grade ensures quality and repeatability every time you use it. The light blue, indigo, and magenta dyes will clearly show the migration of each run, and the dyes migrate independently from the samples allowing the estimation of nucleic acids.

PHYSICAL SPECIFICATIONS

CAS#: Blend
Abs. @ 525nm (1:500 water): 0.3 - 0.6
Abs. @ 588-594nm (1:500 water): 0.5 - 0.8
Abs. @ 635-641nm (1:500 water): 0.5 - 0.7

MOLECULAR BIOLOGY SPECIFICATIONS

DNase assay: None Detected
RNase assay: None Detected
Protease assay: None Detected

STORAGE

- Store at room temperature, Keep tightly sealed, Protect from moisture.



ORDERING INFORMATION

CATALOG# DESCRIPTION

IB01010 6X LOADING DYE-5ML

APPLICATIONS

- IBI's 6X Loading Dye contains 15% Ficoll in a special Tris dye, it is ideal for DNA or RNA gels. The 6X Loading Dye contains three (3) tracking dyes that make estimating sample migration easy and reliable.
 - Dye #1 - Light blue dye is slightly slower than Xylene Cyanol, migrating at around 4 base pairs in a 1% agarose gel.
 - Dye #2 - Indigo dye migrates similar to Bromophenol Blue, at around 600 base pairs in a 1% agarose gel.
 - Dye #3 - Magenta dye migrates at around 150 base pairs in a 1% agarose gel.

PHYSICAL SPECIFICATIONS

CAS#: Blend

MOLECULAR BIOLOGY SPECIFICATIONS

DNase assay: None Detected
RNase assay: None Detected
Protease assay: None Detected

STORAGE

- Store at +4°C, Keep tightly sealed, Protect from moisture.

ORDERING INFORMATION

CATALOG# DESCRIPTION

IB01015 5X RNA GEL LOADING KIT-100 REACTIONS

APPLICATIONS

- IBI's 5X RNA Gel Loading Kit contains reagents for denaturing & loading RNA samples onto a formaldehyde gel, using MOPS as a buffer. The reagents are nuclease free and contain ultra pure deionized formamide & DEPC water.
- The RNA sample is dissolved in 10µl of DEPC water and mixed with 35µl of denaturing solution. Heat the sample to 65°C for 5 minutes. Once the solution has cooled, add 5µl of loading dye. The sample is now ready to load into the gel.

PHYSICAL SPECIFICATIONS

CAS#: Blend
Abs. @ 525nm (1:250 water): 0.15 - 0.35
Abs. @ 588-594nm (1:250 water): 0.32 - 0.47
Abs. @ 635-641nm (1:250 water): 0.25 - 0.40
Contains 10% Glycerol / 2% SDS

MOLECULAR BIOLOGY SPECIFICATIONS

DNase assay: None Detected
RNase assay: None Detected
Protease assay: None Detected

STORAGE

- Store at room temperature, Keep tightly sealed, Protect from moisture.

ORDERING INFORMATION

CATALOG# DESCRIPTION

IB01190 2X PROTEIN LOADING DYE-5ML

APPLICATIONS

- IBI's 2X Protein Loading Dye helps track the migration progression of your sample during polyacrylamide electrophoresis. The loading dye not only gives color to the sample, but also increases the density to ensure efficient distribution in each well. IBI's loading dye migrates independently of the samples, making it easier to estimate the migration of proteins.
- Use 5-10µl per well.

PHYSICAL SPECIFICATIONS

CAS#: 2650-17-1
Molecular Formula: C₂₅H₂₇N₂O₆S₂Na
Molecular Weight: 538.6
Dye Content: Min. 75.0%
Moisture: <10%
Solubility (1% water): Max. 10%
Em (614nm, methanol): Min. 30,000

MOLECULAR BIOLOGY SPECIFICATIONS

DNase assay: None Detected
RNase assay: None Detected
Protease assay: None Detected

ORDERING INFORMATION

CATALOG# DESCRIPTION

IB72120 XYLENE CYANOL FF-25GM

STORAGE

- Store at room temperature, Keep tightly sealed, Protect from moisture.

APPLICATIONS

- Xylene Cyanol FF is used as a tracking dye to monitor the progress of electrophoresis separations. The tracking dye typically migrates with the DNA molecules around 5kb. Therefore, if you are monitoring the progress of longer electrophoresis runs, Xylene Cyanol FF is the tracking dye of choice!

ORDERING INFORMATION

CATALOG#	DESCRIPTION
IB40060	ETHIDIUM BROMIDE-5GM

APPLICATIONS

- Ethidium Bromide can be used to detect both single and double stranded nucleic acids (both DNA and RNA). However, Ethidium Bromide does not intercollate into single-stranded DNA very well; therefore, the fluorescence of single-stranded DNA is relatively poor in comparison to double-stranded DNA. The recommended final concentration of Ethidium Bromide in an agarose gel and/or electrophoresis buffers is 0.5µg/ml.

PHYSICAL SPECIFICATIONS

CAS#:	1239-45-8
Molecular Formula:	C ₂₁ H ₂₀ BrN ₃
Molecular Weight:	394.32
Purity:	>98.0%
Moisture:	Max. 5.0%
Ammonium Bromide:	Max. 0.1%
Melting Point:	272.5 ± 2.5°C

MOLECULAR BIOLOGY SPECIFICATIONS

DNase assay:	None Detected
RNase assay:	None Detected

STORAGE

- Store at room temperature, Keep tightly sealed, Protect from moisture, Light sensitive.

ORDERING INFORMATION

CATALOG#	DESCRIPTION
IB40075	ETHIDIUM BROMIDE SOLUTION-10ML

1 DROP = 0.25µl



ETHIDIUM BROMIDE SOLUTION

PHYSICAL SPECIFICATIONS

CAS#:	1239-45-8
Purity:	Min. 98.0%
Moisture:	Max. 8.0%
Solution Concentration:	10mg per ml

MOLECULAR BIOLOGY SPECIFICATIONS

DNase assay:	None Detected
RNase assay:	None Detected

STORAGE


- Store at room temperature, Keep tightly sealed, Protect from moisture.

APPLICATIONS

- Ethidium Bromide can be used to detect both single and double stranded nucleic acids (both DNA and RNA). However, Ethidium Bromide does not intercollate into single-stranded DNA very well; therefore, the fluorescence of single-stranded DNA is relatively poor in comparison to double-stranded DNA. The recommended final concentration of Ethidium Bromide in an agarose gel and/or electrophoresis buffers is 0.5µg/ml.

ORDERING INFORMATION

CATALOG#	DESCRIPTION
IB01100	PONCEAU S-50GM



PONCEAU S

PHYSICAL SPECIFICATIONS

CAS#:	6226-79-5
Molecular Formula:	C ₂₂ H ₁₂ N ₄ Na ₄ O ₁₃ S ₄
Molecular Weight:	760.61
Loss on drying (%):	≤6%
Em (520nm, methanol):	≥27,000

MOLECULAR BIOLOGY SPECIFICATIONS

DNase assay:	None Detected
RNase assay:	None Detected
Protease assay:	None Detected

STORAGE

- Store at room temperature, Keep tightly sealed, Protect from moisture.

APPLICATIONS

- Ponceau S is a stain used to visualize proteins on nitrocellulose or PVDF membranes. **To create Ponceau S solution;** dissolve 0.5gm of Ponceau S in 1ml of glacial acetic acid, and bring to 100ml with ddi water. **To stain with Ponceau S;** remove blot and rinse with ddi water, stain with Ponceau S solution for 1 to 2 minutes. Rinse for 3 to 5 minutes with ddi water. Red bands should appear against the white background. To reversal stain with Ponceau S Solution; place the nitrocellulose membrane in Ponceau S solution for 5 minutes at room temperature, destain in ddi water for 2 minutes and mark the molecular weight standard with indelible ink, finally to complete the destain soak in water for 10 minutes. Alternatively, use a 0.1N NaOH solution for 5 minutes if proteins of interest are not basic.

PHYSICAL SPECIFICATIONS

CAS#: Blend
 Appearance: Amber Solution
 pH @ 25°C: 0.4 - 0.8
 Abs. @ 465nm: 0.7 - 1.0

PROCEDURE

- 15min. Pre-Wash with Deionized Water (for denaturing gel)
- 60min. Staining with Blue-Clean Protein (completely submerged the gel)
- 30min. Safer Alternative, No more Methanol or Acetic Acid.

STORAGE




- Store at room temperature, Keep tightly sealed, Protect from moisture, Light sensitive.

APPLICATIONS

- Hands off simple staining/destaining procedure
- High sensitivity. Below 20ng of protein per band
- Safer alternative, no methanol or acetic acid required

PROTOCOL

- Fix gels for 2 hours or overnight in the following:
 - 50ml Methanol
 - 12ml Acetic Acid
 - 37% Formaldehyde
 - Glacial Acetic Acid
 - Distilled/Deionized water
- Wash three times for 20 minutes in 35% Ethanol
- Add 10ml of sensitizer solution to 90ml of distilled or deionized water (1:10) and incubate for 2 minutes.
- Wash three times for 5 minutes in distilled or deionized water.
- Add 10ml of Silver Stain Solution to 90ml of water (1:10) and incubate for 20 minutes.
- Wash two times for one minute in distilled or deionized water.
- Add 20ml of Developer Solution to 80ml if distilled or deionized water (1:5) and incubate until bands appear (approximately 10 minutes).
- Stop reaction for 20 minutes in the following solution:
 - 50ml Methanol
 - 12ml Acetic Acid
 - 38ml distilled or deionized water
- Store gels in 1% acetic acid solution at +4°C.

	Snowflake	- Frozen Storage Required
	Thermometer	- Refrigerated Storage Required
	Diamond	- Hazardous Materials

ORDERING INFORMATION

CATALOG# DESCRIPTION

IB01034 BLUE-CLEAN PROTEIN STAIN-1L

Coomassie Blue Staining of proteins in SDS-PAGE gels is a daily procedure in many laboratories. It is very popular among life science researchers due to its good sensitivity and relative ease of use. Traditionally, coomassie blue staining requires a methanol and acetic acid solution to achieve staining and destaining. This process increases the risk for hazardous exposure and produces an unpleasant pungent odor. Blue-Clean Stain is a convenient alternative to traditional coomassie blue staining procedures. Environmentally friendly, this ready-to-use stain does not contain methanol or acetic acid, and does not require hazardous solvents for destaining. This simple "hands-off" staining/destaining procedure saves valuable time while reducing the handling of hazardous materials and solvent waste in your lab. Blue-Clean staining exhibits sensitivity below 20ng of protein per band. Packaged as a 1X ready-to-use solution, Blue-Clean Stain only requires water for the prewashing and destaining procedures.



ORDERING INFORMATION

CATALOG# DESCRIPTION

IB01037 SILVER STAIN KIT

Techniques involving staining of nucleic acids and proteins using ethidium bromide and coomassie blue stain are useful for a wide spectrum of applications, but both methods are limited in their sensitivity. When very low concentrations of nucleic acids or proteins need to be stained, silver staining methods still provide the most effective alternative. IBI's Silver Stain Kit is a rapid and easy-to-use method to achieve excellent silver staining of nucleic acids or proteins in polyacrylamide gels. This method requires only one hour after fixing the gel. Silver Stain has outstanding sensitivity with very low background. IBI's Silver Stain Kit also provides sufficient, pre-measured reagents for staining 20 gels, along with complete easy to follow directions for use. All reagents are room temperature stable for more than one year, reconstituted silver binding agent should be refrigerated for long term storage.



ORDERING INFORMATION

CATALOG#	DESCRIPTION
IB74020	ACRIDINE ORANGE-25GM

APPLICATIONS

- Acridine Orange is a metachromatic stain that can be used to detect both single and double stranded nucleic acids. Acridine orange interacts with polynucleotides, either by intercalation between the stacked bases of dsDNA or RNA, or when the polynucleotide is predominantly single stranded by binding electrostatically to the phosphate backbone to produce a stacked array. If the stain intercalates between the stacked bases of the double stranded nucleic acid it will fluoresce green at 530nm. If acridine orange binds to the single stranded nucleic acid, it will fluoresce orange-red at 640nm. Therefore, acridine orange is a very informative reagent when denaturing the structure of the nucleic acids that have undergone electrophoresis. Acridine orange is also more sensitive than ethidium bromide if the electrophoresis system has been denatured by glyoxalation. For a 1.5% agarose gel, stain with 30mg/ml of acridine orange in 10mM sodium phosphate (pH 7.0) for 15 minutes at room temperature in the dark. To destain, run hot tap water over the gel for 15-20 minutes.

ORDERING INFORMATION

CATALOG#	DESCRIPTION
IB74040	BROMOPHENOL BLUE-25GM

APPLICATIONS

- Bromophenol Blue is used as a tracking dye to monitor the progress of electrophoresis separations. Bromophenol Blue usually comigrates with DNA molecules around 0.5kb. Therefore, bromophenol blue provides an index of the mobility of the fastest fragments and is valuable in determining the length of the gel, over which the separation of DNA has occurred.

ORDERING INFORMATION

CATALOG#	DESCRIPTION
IB74050	METHYLENE BLUE-25GM

APPLICATIONS

- Methylene Blue is a recommended stain for agarose gels when recovery of DNA is required. The benefit of using methylene blue (rather than ethidium bromide) is the ability to use a white light box for visualization and photography in place of a UV light box. Using a white light box eliminates mutations that may affect the integrity of the isolated DNA.



PHYSICAL SPECIFICATIONS

CAS#:	10127-02-3
Molecular Formula:	$C_{17}H_{14}ClN_3 \cdot \frac{1}{2}ZnCl_2$
Molecular Weight:	370.0
Em (490nm water):	Min. 36,500

STORAGE

- Store at room temperature, Keep tightly sealed, Protect from moisture, Light sensitive.

PHYSICAL SPECIFICATIONS

CAS#:	115-39-9
Molecular Formula:	$C_{19}H_{16}Br_4O_2S$
Molecular Weight:	669.96
Purity:	>98.0%
pH (1.0M water):	7.0 \pm 0.4
Insolubles:	<0.01%

STORAGE

- Store at room temperature, Keep tightly sealed, Protect from moisture, Light sensitive.

PHYSICAL SPECIFICATIONS

CAS#:	7220-79-3
Molecular Formula:	$C_{16}H_{18}ClN_3S \cdot 3H_2O$
Molecular Weight:	373.9
Em (661nm water):	>69,000
Arsenic:	0.0008%
Copper:	<0.01%
Moisture:	Max. 15.0%
Zinc:	<0.001%
Lead:	<0.001%
Dye Content (Dry Basis):	98 - 103%

STORAGE

- Store at room temperature, Protect from moisture.

STAINING DNA PROTOCOL

- 1.) Stain a 1% gel for 15 minutes with 0.02% Methylene Blue Solution in distilled water.
- 2.) Destain with distilled water for 15 minutes.

STAINING RNA PROTOCOL (bonded to nylon or nitrocellulose membranes)

- 1.) Soak the dried filter in 5% acetic acid at room temperature for 15 minutes.
- 2.) Place the filter in 0.5M sodium acetate (pH 5.2) and 0.04% methylene blue for 5 - 10 minutes at room temperature.
- 3.) Rinse the filter in water for 5 - 10 minutes.

STORAGE

- Store at +4°C, Keep tightly sealed, Protect from moisture, Light sensitive.

APPLICATIONS

- Staining proteins in analytical work
- Staining 2-D gels
- Staining gels for long-term documentation

STAINING PROCEDURE

- Remove stacking gel
- 10min. Rinse gel with 300ml of water(for denaturing gels) three (3) times
- 60min. Staining with See Band (completely submerge gel) 20ml-mini gel
- 30min. Destain in deionized water (Do not continue overnight)

NOTE: Longer destaining may be necessary for optimal results.

TIME TO REMOVE SDS

- <14% 10min. in 300ml water, three (3) times
- 14%-16% 30min. in 300ml water, three (3) times
- >16% 60min. in 300ml water, three (3) times

RECOMMENDED USE

- Clear and sharp bands
- High sensitivity below 10ng of protein per band
- Track staining intensity during staining
- Low background
- Washes with only water
- Stain does not interfere with Western Blottings
- Storage of 12 months @ +4°C

Seeband Protein Stain is a convenient alternative to traditional coomassie blue staining procedures. *Environmentally friendly*, this ready-to-use stain does not contain methanol or acetic acid, and does not require hazardous solvents for destaining. This simple “hands-off” staining/destaining procedure saves valuable time while reducing the handling of hazardous materials and solvent waste in your laboratory.

ORDERING INFORMATION

CATALOG#	DESCRIPTION
IB01033	SEEBAND PROTEIN STAIN (ANALYTICAL GRADE)-500ML

STORAGE

- Store at +4°C, Keep tightly sealed, Protect from moisture, Light sensitive.

APPLICATIONS

- Staining proteins in extraction work
- Staining gel prior to Western Blotting
- Staining 2-D gels
- Staining gels for MALDI-MS

STAINING PROCEDURE

- Remove stacking gel
- 10min. Rinse gel with 300ml of water(for denaturing gels) three (3) times
- 60min. Staining with See Band (completely submerge gel) 20ml-mini gel
- 30min. Destain in deionized water (Do not continue overnight)

NOTE: Longer destaining may be necessary for optimal results.

TIME TO REMOVE SDS

- <14% 10min. in 300ml water, three (3) times
- 14%-16% 30min. in 300ml water, three (3) times
- >16% 60min. in 300ml water, three (3) times




RECOMMENDED USE

- Clear and sharp bands
- High sensitivity below 10ng of protein per band
- Track staining intensity during staining
- Low background
- Washes with only water
- Stain does not interfere with Western Blottings
- Storage of 12 months @ +4°C

Seeband Protein Stain is a convenient alternative to traditional coomassie blue staining procedures. *Environmentally friendly*, this ready-to-use stain does not contain methanol or acetic acid, and does not require hazardous solvents for destaining. This simple “hands-off” staining/destaining procedure saves valuable time while reducing the handling of hazardous materials and solvent waste in your laboratory.

ORDERING INFORMATION

CATALOG#	DESCRIPTION
IB01035	SEEBAND PROTEIN STAIN (EXTRACTION GRADE)-500ML

	Snowflake	-	Frozen Storage Required
	Thermometer	-	Refrigerated Storage Required
	Diamond	-	Hazardous Materials

ENZYMES / SUBSTRATES

REAGENTS

PHYSICAL SPECIFICATIONS

CAS#:	3483-12-3
Molecular Formula:	C ₄ H ₁₀ O ₂ S ₂
Molecular Weight:	154.25
-SH Content:	Min. 99.4%
A ₂₈₃ (0.02M water):	Max. 0.05
-S-S Content:	Max. 0.5%
Melting Point:	41 ± 2°C
IR:	Pass
Moisture:	Max. 0.5%
Solubility (5% water):	Clear, Colorless

MOLECULAR BIOLOGY SPECIFICATIONS

DNase assay:	None Detected
RNase assay:	None Detected
Protease assay:	None Detected

STORAGE

- Store at -20°C, Keep tightly sealed, Warm to room temperature before opening. A DTT solution should be stored in small aliquots at -20°C. It is advisable that DTT solutions should be made fresh prior to use.
- Note:** Do NOT autoclave DTT or solutions containing DTT.

ORDERING INFORMATION

CATALOG# DESCRIPTION

IB21040	DTT (DITHIOTHREITOL-CLELAND'S REAGENT)-5GM
IB21045	DTT (DITHIOTHREITOL-CLELAND'S REAGENT)-25GM

APPLICATIONS

- DTT is used in enzyme preparations to prevent formation of disulfides from sulfhydryl groups, thereby stabilizing the enzyme.



PHYSICAL SPECIFICATIONS

CAS#:	56-81-5
Molecular Formula:	C ₃ H ₈ O ₃
Molecular Weight:	92.09
Purity:	Min. 99.5%
Sulfate:	<0.0001%
Arsenic:	<0.0003%
Chlorine:	<0.0003%
Iron:	<0.0005%
Residue on ignition:	≤0.01
Density:	1.249 - 1.262gm/ml
Specific Gravity:	1.2635 - 1.2615

MOLECULAR BIOLOGY SPECIFICATIONS

DNase assay:	None Detected
RNase assay:	None Detected
Protease assay:	None Detected

APPLICATIONS

- Glycerol is used to prevent freezing during enzyme and bacterial preparations.

STORAGE

- Store at room temperature, Protect from moisture, Keep tightly sealed.

ORDERING INFORMATION

CATALOG# DESCRIPTION

IB15760	GLYCEROL-500ML
IB15762	GLYCEROL-1L



PHYSICAL SPECIFICATIONS

CAS#:	39450-01-6
Molecular Weight:	~18,000 Da
Activity (@37°C U/mg):	Min. 30U/mg

MOLECULAR BIOLOGY SPECIFICATIONS

DNase assay:	None Detected
RNase assay:	None Detected

STORAGE

- Store at -20°C, Warm to room temperature before opening, Keep tightly sealed.


APPLICATIONS

- Proteinase K is a stable >50mg/ml (powder) highly active proteolytic enzyme that is purified from the mold tritirachium album. The enzyme has two binding sites for Ca⁺⁺, which lie some distance from the active site, and is not directly involved in the catalytic mechanism. The high activity of the enzyme solution (0.6U/μl) means that you can expect more complete digestion of your sample.

ORDERING INFORMATION

CATALOG# DESCRIPTION

IB05400	PROTEINASE K-100MG
---------	--------------------



TO MAKE A 20MG/ML PROTEINASE K SOLUTION (CONC.)

Tris	10mM
Calcium Chloride	1mM
Glycerol	30%
pH adjusted to 8.0 w/HCl	

ENZYMES / SUBSTRATES

REAGENTS

ORDERING INFORMATION

CATALOG# DESCRIPTION

IB02100	IPTG-1GM
IB02105	IPTG-5GM
IB02125	IPTG-25GM

PHYSICAL SPECIFICATIONS

CAS#:	367-93-1
Molecular Formula:	C ₉ H ₁₈ O ₅ S
Molecular Weight:	238.31
Purity:	Min. 99.0%
Melting Point:	110 - 114°C
IR (KBr):	Conforming as Standard
Dioxane (GLC):	<1ppm

MOLECULAR BIOLOGY SPECIFICATIONS

Colorimetric Transformation assay: Pass

STORAGE

- Store desiccated at -20°C, Keep tightly sealed, Warm to room temperature before opening.

APPLICATIONS

- IPTG is an inducer of β -galactosidase activity in *E. coli*. It is commonly used in conjunction with X-GAL to detect lac gene expression in cloning applications, thus allowing detection of recombinant molecules.

PROTOCOL FOR STOCK SOLUTION

- Prepare a stock solution of IPTG at 0.1M (23.8mg/ml) in sterile distilled water and filter sterilize through a 0.22 micron disposable filter. Add 3ml per liter of IPTG stock solution and any of the solutions listed to the left. When pouring plates add 15gm bacto-agar per liter.

10X M9 Salts Solution:

60gm	Ma ₂ HPO ₄
30gm	KH ₂ PO ₄
5gm	NaCl
10gm	NH ₄ Cl

Adjust the pH to 7.4, add water to a final volume of 1L, & autoclave.

For M9 media add:

100ml	10X M9 Salt Solution
1ml	0.01M CaCl ₂
2ml	1.0M MgSO ₄

Adjust to a final volume of 1L with distilled water.

Supplement with (optional):

0.2%	Casamino Acids
0.4%	Glucose
40mg/ml	Thiamine HCl

2X YT Media:

16gm/L	Tryptone
10gm/L	Yeast Extract
5gm/L	NaCl



IPTG

ORDERING INFORMATION

CATALOG# DESCRIPTION

IB01090	PMSF (PHENYLMETHYL SULFONYL FLUORIDE)-5GM
---------	---

APPLICATIONS

- An irreversible serine protease inhibitor of trypsin, chymotrypsin, kallikrein, subtilisin, and thrombin, as well as cysteine protease papain.

Working Range

0.1 - 1.0mM
0.2 Soluble in Ethanol and Methanol

PHYSICAL SPECIFICATIONS

CAS#:	329-98-6
Molecular Formula:	C ₈ H ₇ CH ₂ SO ₂ F
Molecular Weight:	174.02
Purity:	Min. 99.0%
Melting Point:	91 - 95°C
Solubility(5% EtOH):	Pass
Identification:	Pass
TLC:	One spot

MOLECULAR BIOLOGY SPECIFICATIONS

DNase assay:	None Detected
RNase assay:	None Detected
Protease assay:	None Detected

STORAGE

- Store at room temperature, Keep tightly sealed, Protect from moisture.

ORDERING INFORMATION

CATALOG# DESCRIPTION

IB05406	PROTEINASE K SOLUTION-5ML
---------	---------------------------

STORAGE

- Store at +4°C, Keep tightly sealed.



PHYSICAL SPECIFICATIONS

CAS#:	Blend
Solution Components	
Tris:	10mM
Calcium Chloride:	1mM
Glycerol:	30%
pH adjusted to 8.0 with HCl	

MOLECULAR BIOLOGY SPECIFICATIONS

DNase assay (Endo & Exo):	None Detected
RNase assay:	None Detected
Proteinase activity on gel substrate:	Pass

APPLICATIONS

- IBI's Proteinase K Solution is a stable 20mg/ml highly active proteolytic enzyme that is purified from the mold *Tritirachium album*. The enzyme has two binding sites for Ca ++, which lie some distance from the active site and is not directly involved in the catalytic mechanism. The high activity of the enzyme solution (0.6U/μl) means that you can keep the reaction volumes low, as low as 200mg/ml. It also means that you can expect more complete digestion of your samples. Examples: For a final concentration of 200μg/ml, add 1ml of Proteinase K Solution to your sample.

ENZYMES / SUBSTRATES

REAGENTS

PHYSICAL SPECIFICATIONS

CAS#: 7240-90-6
Molecular Formula: $C_{14}H_{15}BrClNO_6$
Molecular Weight: 408.64
Purity: Min. 99%
Solubility (2% in N,N' DMF): Clear, Colorless

MOLECULAR BIOLOGY SPECIFICATIONS

Optical rotation (C=1): Pass

STORAGE

- Store at -20°C, Keep tightly sealed, Warm to room temperature prior to opening.

APPLICATIONS

- X-GAL is histochemical substrate for β -galactosidase, allowing the detection of lac gene expression. To make a stock solution of X-GAL, dissolve the X-GAL in dimethylformamide to make a 20mg/ml solution. **Note:** Use either a glass or polypropylene tube. Add 4ml of this solution per liter of culture media. The tube containing the solution should be wrapped in aluminum foil to prevent damage from light, and should be stored at -20°C. **Note:** It is not necessary to sterilize X-GAL solutions by filtration.

ORDERING INFORMATION



CATALOG# DESCRIPTION

IB02260 X-GAL-1GM
IB02264 X-GAL-100MG



PHYSICAL SPECIFICATIONS

CAS#: Blend

MOLECULAR BIOLOGY SPECIFICATIONS

Name	MW	Description
AEBSF	239.5	Irreversible serine protease inhibitor
Aprotinin	6512.0	Competitive reversible serine protease inhibitor
Bestatin	308.4	Competitive inhibitor of aminopeptidases
E64	357.4	Irreversible inhibitor cysteine proteases
Leupeptin	493.6	Reversible inhibitor of cysteine proteases

STORAGE

- Store at room temperature, Keep tightly sealed, Protect from moisture, Freeze when reconstituted.

ORDERING INFORMATION

CATALOG# DESCRIPTION

IB01070 PROTEASE INHIBITOR COCKTAIL-1ML

APPLICATIONS

- Protease Inhibitor Cocktail is an easy to use lyophilized power. Add 1ml of ddH₂O to the vial and mix gently. The reconstituted protease cocktail should be aliquoted into multiple tubes and frozen. This material can be stored up to two weeks. A 1ml vial of protease cocktail is sufficient for 100ml of sample buffer. See components of protease cocktail listed below:

Name	MW	Amt/1ml	Working Dilution
AEBSF	239.5	50mM	0.5mM
Aprotinin	6512.0	30 μ M	0.3 μ M
Bestatin	308.4	1mM	10.0 μ M
E64	357.4	1mM	10.0 μ M
Leupeptin	493.6	1mM	10.0 μ M

PHYSICAL SPECIFICATIONS

CAS#: Blend

MOLECULAR BIOLOGY SPECIFICATIONS

Name	MW	Description
AEBSF	239.5	Irreversible serine protease inhibitor
Aprotinin	6512.0	Competitive reversible serine protease inhibitor
Bestatin	308.4	Competitive inhibitor of aminopeptidases
E64	357.4	Irreversible inhibitor cysteine proteases
Leupeptin	493.6	Reversible inhibitor of cysteine proteases
EDTA	372.24	Reversible metalloprotease inhibitor

STORAGE

- Store at room temperature, Keep tightly sealed, Protect from moisture, Freeze when reconstituted.

ORDERING INFORMATION

CATALOG# DESCRIPTION

IB01150 PROTEASE INHIBITOR COCKTAIL W/EDTA-1ML

APPLICATIONS

- Protease Inhibitor Cocktail w/EDTA is an easy to use lyophilized power. Add 1ml of ddH₂O to the vial and mix gently. The reconstituted protease cocktail w/EDTA should be aliquoted into multiple tubes and frozen. This material can be stored up to two weeks. A 1ml vial of protease cocktail is sufficient for 100ml of sample buffer. See components of protease cocktail listed below:

Name	MW	Amt/1ml	Working Dilution
AEBSF	239.5	50mM	0.5mM
Aprotinin	6512.0	30 μ M	0.3 μ M
Bestatin	308.4	1mM	10.0 μ M
E64	357.4	1mM	10.0 μ M
Leupeptin	493.6	1mM	10.0 μ M
EDTA	372.24	1mM	1.0-10.0 μ M

	Snowflake	- Frozen Storage Required
	Thermometer	- Refrigerated Storage Required
	Diamond	- Hazardous Materials