

IB50000 - QSH LAB-PAL COMBS

CATALOG#	DESCRIPTION	WELL WIDTH	SAMPLE VOLUME PER MM GEL
IB50040	1.0mm x 5 tooth	6.0mm	6.0µl
IB50050	1.0mm x 10 tooth	2.6mm	2.6µl
IB50080	1.5mm x 5 tooth	6.0mm	9.0µl
IB50090	1.5mm x 10 tooth	2.6mm	3.9µl
IB50060	1.0mm, 0 marker, 1 sample	44.0mm	44.0µl
IB50070	1.0mm, 1 marker, 1 sample	6.0, 34.6mm	6.0, 34.6µl
IB50100	1.5mm, 0 marker, 1 sample	44.0mm	66.0µl
IB50110	1.5mm, 1 marker, 1 samples	6.0, 34.6mm	9.0, 52.0µl

IB51000 - QS-710 COMBS

CATALOG#	DESCRIPTION	WELL WIDTH	SAMPLE VOLUME PER MM GEL
IB51040	1.0mm x 8 tooth	4.7mm	4.7µl
IB51045	1.0mm x 12 tooth	2.8mm	2.8µl
IB51050	1.0mm x 15 tooth	2.2mm	2.2µl
IB51055	1.5mm x 8 tooth	4.7mm	7.0µl
IB51060	1.5mm x 12 tooth	2.8mm	4.0µl
IB51065	1.5mm x 15 tooth	2.2mm	3.3µl
IB51046	2.0mm x 12 tooth	2.8mm	6.0µl
IB51090	2.0mm x 10 tooth	4.0mm	8.0µl
IB51095	3.0mm x 6 tooth	5.8mm	17.3µl
IB51070	1.0mm, 0 marker, 1 sample	62.6mm	62.6µl
IB51075	1.0mm, 1 marker, 1 sample	4.7, 54.7mm	4.7, 54.7µl
IB51080	1.5mm, 0 marker, 1 sample	62.6mm	94.0µl
IB51085	1.5mm, 1 marker, 1 sample	4.7, 54.7mm	7.0, 82.0µl

IB53000 / IB52000 - MP-1015 / MPH COMBS

CATALOG#	DESCRIPTION	WELL WIDTH	SAMPLE VOLUME PER MM GEL
IB52060	1.0mm x 10 tooth	6.3mm	6.3µl
IB52090	1.0mm x 16 tooth	3.5mm	3.5µl
IB52160	1.0mm x 20 tooth	2.6mm	2.6µl
IB52165	1.5mm x 20 tooth	2.6mm	4.5µl
IB52155	1.5mm x 22 tooth	2.6mm	3.9µl
IB52070	2.0mm x 10 tooth	6.3mm	12.6µl
8920316	2.0mm x 11 tooth, MCP**	6.3mm	12.6µl
IB53070	2.0mm x 13 tooth	4.0mm	8.0µl
IB52100	2.0mm x 16 tooth	3.5mm	7.0µl
IB52170	2.0mm x 20 tooth	2.6mm	5.2µl
IB52080	3.0mm x 10 tooth	6.3mm	18.9µl
8695108	3.0mm x 11 tooth, MCP**	6.3mm	18.9µl
IB52110	3.0mm x 16 tooth	3.5mm	10.5µl
IB52180	3.0mm x 20 tooth	2.6mm	7.8µl
IB52120	1.0mm, 2 markers, 1 sample	2.9, 81.8mm	2.9, 82.0µl
IB52130	2.0mm, 2 markers, 1 sample	2.9, 81.8mm	5.8, 163.0µl
IB53080	3.0mm, 2 markers, 1 sample	2.9, 81.8mm	8.7, 245.0µl
IB52140	1.0mm, 2 markers, 2 samples	2.9, 39.3mm	2.9, 39.0µl
IB52150	2.0mm, 2 markers, 2 samples	2.9, 39.3mm	5.8, 78.0µl
IB53090	3.0mm, 2 markers, 2 samples	2.9, 39.3mm	8.7, 117.0µl

IB56000 / IB54000 - HR-2025 / HRH COMBS

CATALOG#	DESCRIPTION	WELL WIDTH	SAMPLE VOLUME PER MM GEL
IB54180	1.0mm x 10 tooth	16.6mm	16.6µl
IB56050	1.0mm x 12 tooth	13.1mm	13.1µl
IB56055	1.0mm x 15 tooth	9.8mm	9.8µl
IB54060	1.0mm x 20 tooth	6.5mm	6.5µl
IB54090	1.0mm x 30 tooth	4.0mm	4.0µl
IB56100	1.0mm x 40 tooth	3.2mm	3.2µl
IB56101	1.0mm x 44 tooth, MCP**	2.9mm	2.9µl
IB54190	2.0mm x 10 tooth	16.6mm	49.8µl
IB56115	2.0mm x 12 tooth	13.1mm	26.2µl
IB56125	2.0mm x 15 tooth	9.8mm	19.6µl
IB54070	2.0mm x 20 tooth	6.5mm	13.0µl
8091100	2.0mm x 22 tooth, MCP**	6.3mm	12.6µl

MCP** - Multichannel Pipette Compatible

IB56000 / IB54000 - HR-2025 / HRH COMBS

CATALOG#	DESCRIPTION	WELL WIDTH	SAMPLE VOLUME PER MM GEL
IB54100	2.0mm x 30 tooth	4.0mm	8.0µl
IB56135	2.0mm x 40 tooth	3.2mm	6.4µl
IB56102	2.0mm x 44 tooth, MCP**	2.9mm	5.8µl
IB54080	3.0mm x 20 tooth	6.5mm	19.5µl
8093767	3.0mm x 22 tooth, MCP**	6.3mm	18.9µl
IB54110	3.0mm x 30 tooth	4.0mm	12.0µl
IB56140	3.0mm x 40 tooth	3.2mm	9.6µl
IB56103	3.0mm x 44 tooth, MCP**	2.9mm	8.7µl
IB54120	1.0mm, 2 markers, 1 sample	4.0, 184.0mm	4.0, 184.0µl
IB54140	1.0mm, 2 markers, 2 samples	4.0, 91.0mm	4.0, 91.0µl
IB54130	2.0mm, 2 markers, 1 sample	4.0, 184.0mm	8.0, 368.0µl
IB54150	2.0mm, 2 markers, 2 samples	4.0, 91.0mm	8.0, 182.0µl
IB56145	3.0mm, 2 markers, 1 sample	4.0, 184.0mm	12.0, 552.0µl
IB56150	3.0mm, 2 markers, 2 samples	4.0, 91.0mm	12.0, 273.0µl

IB57000 - HR2525 COMBS

CATALOG#	DESCRIPTION	WELL WIDTH	SAMPLE VOLUME PER MM GEL
IB57125	1.0mm x 25 tooth, MCP**	5.8mm	5.8µl
IB57150	1.0mm x 50 tooth, MCP**	2.9mm	2.9µl
IB57225	2.0mm x 25 tooth, MCP**	5.8mm	11.6µl
IB57250	2.0mm x 50 tooth, MCP**	2.9mm	5.8µl
IB57325	3.0mm x 25 tooth, MCP**	5.8mm	17.5µl
IB57350	3.0mm x 50 tooth, MCP**	2.9mm	8.7µl
IB57425	1.5mm x 25 tooth, MCP**	5.8mm	8.7µl
IB57450	1.5mm x 50 tooth, MCP**	2.9mm	4.4µl

IB80000 - STS-45i COMBS

CATALOG#	DESCRIPTION	WELL WIDTH	SAMPLE VOLUME PER MM GEL
IB80040	0.2mm x 32 tooth	8.0mm	1.6µl
IB80050	0.2mm x 64 tooth	3.0mm	0.6µl
IB80100	0.4mm x 32 tooth	8.0mm	3.2µl
IB80110	0.4mm x 64 tooth	3.0mm	1.2µl
IB80250	0.2mm x 64 sharkstooth	3.0mm	0.6µl
IB80355	0.2mm x 64&96 sharkstooth	3.0, 2.0mm	0.2, 0.4µl
IB80310	0.4mm x 64 sharkstooth	3.0mm	1.2µl
IB80350	0.2mm x 96 sharkstooth	2.0mm	0.4µl
IB80380	0.4mm x 96 sharkstooth	2.0mm	0.8µl
IB80385	0.4mm x 64&96 sharkstooth	3.0, 2.0mm	1.2, 0.8µl

IB62000 - VCV COMBS

CATALOG#	DESCRIPTION	WELL WIDTH	SAMPLE VOLUME PER MM GEL
IB62045	0.8mm x 12 tooth	10.0mm	40µl
IB62105	0.8mm x 20 tooth	5.0mm	25µl
IB62055	1.5mm x 12 tooth	10.0mm	70µl
IB62115	1.5mm x 20 tooth	5.0mm	75µl
IB62065	3.0mm x 12 tooth	10.0mm	120µl
IB62125	3.0mm x 20 tooth	5.6mm	150µl
IB62345	0.8mm, 2 markers, 2 sample	6.3, 72.0mm	5.0, 57.6µl
IB62355	1.5mm, 2 markers, 2 sample	6.3, 72.0mm	9.5, 108.0µl
IB62365	3.0mm, 2 markers, 1 sample	6.3, 147.0mm	18.9, 441.0µl
IB62375	3.0mm, 2 markers, 2 sample	6.3, 72.0mm	18.9, 216.0µl

IB94000 - MAGELIN COMBS

CATALOG#	DESCRIPTION	WELL WIDTH	SAMPLE VOLUME PER MM GEL
IB92030	0.8mm x 10 tooth	5.9mm	28.1µl
IB92031	0.8mm x 12 tooth	4.2mm	20.3µl
IB92032	0.8mm x 15 tooth	3.4mm	16.2µl
IB92033	0.8mm, 0 markers, 1 sample	7.0mm	336.0µl

JORDAN COMBS

COMBS

JSB-30 - COMBS

CATALOG#	DESCRIPTION	WELL WIDTH	SAMPLE VOLUME PER MM GEL
CSM-0808	0.8mm x 8 tooth	5.8mm	4.7µl
CSM-0813	0.8mm x 13 tooth	4.0mm	3.2µl
CSM-0818	0.8mm x 18 tooth	2.3mm	1.9µl
CSM-1508	1.5mm x 8 tooth	5.8mm	8.7µl
CSM-1513	1.5mm x 13 tooth	4.0mm	6.0µl
CSM-1518	1.5mm x 18 tooth	2.3mm	3.5µl

JSB-302 - COMBS

CATALOG#	DESCRIPTION	WELL WIDTH	SAMPLE VOLUME PER MM GEL
CSD-0816	0.8mm x 16 tooth	5.8mm	4.7µl
CSD-0827	0.8mm x 27 tooth	4.0mm	3.2µl
CSD-0836	0.8mm x 36 tooth	2.3mm	1.9µl
CSD-1516	1.5mm x 16 tooth	5.8mm	8.7µl
CSD-1522	1.5mm x 22 tooth	5.0mm	7.5µl
CSD-1527	1.5mm x 27 tooth	4.0mm	6.0µl
CSD-1533	1.5mm x 33 tooth	3.0mm	4.5µl
CSD-1536	1.5mm x 36 tooth	2.3mm	3.5µl

JSB-96 - COMBS

CATALOG#	DESCRIPTION	WELL WIDTH	SAMPLE VOLUME PER MM GEL
CB96-0825	0.8mm x 25 tooth, MCP**	5.8mm	4.7µl
CB96-0840	0.8mm x 40 tooth	4.0mm	3.2µl
CB96-0850	0.8mm x 50 tooth, MCP**	2.9mm	2.3µl
CB96-1525	1.5mm x 25 tooth, MCP**	5.8mm	8.7µl
CB96-1540	1.5mm x 40 tooth	4.0mm	6.0µl
CB96-1550	1.5mm x 50 tooth, MCP**	2.9mm	4.4µl
CB96-3025	3.0mm x 25 tooth, MCP**	5.8mm	17.4µl
CB96-3040	3.0mm x 40 tooth	4.0mm	12.0µl
CB96-15P	1.5mm, 2 markers, 1 sample	6.4, 204.0mm	9.6, 306.0µl
CB96-30P	3.0mm, 2 markers, 1 sample	6.4, 204.0mm	19.2, 612.0µl

JSB-75 - COMBS

CATALOG#	DESCRIPTION	WELL WIDTH	SAMPLE VOLUME PER MM GEL
CB75-0808	0.8mm x 8 tooth	9.9mm	7.9µl
CB75-0813	0.8mm x 13 tooth	5.8mm	4.7µl
CB75-0822	0.8mm x 22 tooth	4.0mm	16.0µl
CB75-0830	0.8mm x 30 tooth	2.5mm	2.0µl
CB75-1508	1.5mm x 8 tooth	9.9mm	14.9µl
CB75-1513	1.5mm x 13 tooth	5.8mm	8.7µl
CB75-1522	1.5mm x 22 tooth	4.0mm	6.0µl
CB75-1530	1.5mm x 30 tooth	2.5mm	3.7µl
CB75-3008	3.0mm x 8 tooth	9.9mm	29.72µl
CB75-3013	3.0mm x 13 tooth	5.8mm	17.5µl
CB75-3022	3.0mm x 22 tooth	4.0mm	12.0µl

JSB-120 - COMBS

CATALOG#	DESCRIPTION	WELL WIDTH	SAMPLE VOLUME PER MM GEL
CB120-0812	0.8mm x 12 tooth	10.0mm	8.1µl
CB120-0821	0.8mm x 21 tooth	5.8mm	4.7µl
CB120-0836	0.8mm x 36 tooth	4.0mm	3.2µl
CB120-0848	0.8mm x 48 tooth	2.5mm	2.0µl
CB120-1512	1.5mm x 12 tooth	10.0mm	15.2µl
CB120-1521	1.5mm x 21 tooth	5.8mm	8.7µl
CB120-1536	1.5mm x 36 tooth	4.0mm	6.0µl
CB120-1548	1.5mm x 48 tooth	2.5mm	3.7µl
CB120-3012	3.0mm x 12 tooth	10.0mm	30.5µl
CB120-3021	3.0mm x 21 tooth	5.8mm	17.5µl
CB120-3036	3.0mm x 36 tooth	4.0mm	12.0µl

MCP** - Multichannel Pipette Compatible

JVD-80 - COMBS

CATALOG#	DESCRIPTION	WELL WIDTH	SAMPLE VOLUME PER MM GEL
CP-0406	0.4mm x 6 tooth	25.4mm	10.1µl
CP-0412	0.4mm x 12 tooth	9.9mm	4.0µl
CP-0418	0.4mm x 18 tooth	5.8mm	2.3µl
CP-0427	0.4mm x 27 tooth	4.0mm	1.6µl
CP-0806	0.8mm x 6 tooth	25.4mm	20.3µl
CP-0812	0.8mm x 12 tooth	9.9mm	7.9µl
CP-0818	0.8mm x 18 tooth	5.8mm	4.6µl
CP-0827	0.8mm x 27 tooth	4.0mm	3.2µl
CP-1506	1.5mm x 6 tooth	25.4mm	38.1µl
CP-1512	1.5mm x 12 tooth	9.9mm	14.9µl
CP-1518	1.5mm x 18 tooth	5.8mm	8.7µl
CP-1527	1.5mm x 27 tooth	4.0mm	6.0µl
CP-3006	3.0mm x 6 tooth	25.4mm	76.2µl
CP-3012	3.0mm x 12 tooth	9.9mm	29.7µl
CP-3018	3.0mm x 18 tooth	5.8mm	17.4µl
CP-3027	3.0mm x 27 tooth,	4.0mm	12.0µl
CP-04P	0.4mm, 2 markers, 1 sample	6.3, 152.4mm	2.5, 61.0µl
CP-08P	0.8mm, 2 markers, 1 sample	6.3, 152.4mm	5.0, 121.9µl
CP-15P	1.5mm, 2 markers, 1 sample	6.3, 152.4mm	9.5, 228.6µl
CP-30P	3.0mm, 2 markers, 1 sample	6.3, 152.4mm	18.9, 457.2µl
PF-26P	0.4mm x 26 tooth*	3.0mm	1.2µl
PF-36P	0.4mm x 32 tooth*	2.5mm	1.0µl
PF-46P	0.4mm x 46 tooth*	2.0mm	0.8µl

JQT-182 - COMBS

CATALOG#	DESCRIPTION	WELL WIDTH	SAMPLE VOLUME PER MM GEL
CQ-0418	0.4mm x 18 tooth	13.6mm	5.44µl
CQ-0432	0.4mm x 32 tooth	5.8mm	2.32µl
CQ-0454	0.4mm x 54 tooth	4.0mm	1.6µl
PF-64Q	0.4mm x 64 picket fence	1.9mm	0.8µl
PF-80Q	0.4mm x 80 picket fence	1.9mm	0.8µl



IBI ANALYTICAL COMBS



IBI PREPARATORY COMBS